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Please read the instructions for use carefully before using the product.

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1. INDICATION

The Bonder is used as a layer between the alloy and the ceramic. It blocks escaping metal oxides and neutralize differences in the thermal expansion coefficient.

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2. USE

We recommend to test the desired alloy with **Crea Alloy Bond** to ensure successful use.

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2.1. THE METAL FRAMEWORK

Should be prepared according to the manufacturer's instructions.

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2.2. APPLYING

Stir before using! Apply the bonder in a thin layer but mask the framework completely.

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Important: Do not moisten dried Bonder. Close the cap after using.

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2.4. APPEARANCE

After firing the Bonder shows a yellowish, slight shiny surface.

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Note: The colour of the fired Bonder depends on the alloy's composition.

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2.5. OPAQUE FIRING AND FURTHER PROCESSING

1 st Opaque Firing	
Creation CC	NPA 980 °C
Creation LF	NPA 830 °C

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The further processing according the manufacturer's instructions.

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3. SAFETY

- Inhale: Fresh air.
- Skin contact: Wash off using water.
- Eye contact: Rinse eyes under running water with eyelid being opened.
- After swallowing: Drink a lot of water. If complaints persist, seek medical attention.
- Protective equipment: Protective gloves, safety goggles and protective working clothes.

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4. NOTE

The instructions for use are always updated according to further developments. We recommend to read the instructions for use before using a new package.

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5. IMPORTANT

Original packed
By open packing no right to return.

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