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PACKUNGSGRÖSSEN

1. Sortiment

- | | |
|--|--------|
| CREATION VC Make up Neo Kit (Make up Neo stain A, stain B, stain C, white, eggshell, lemon yellow, mandarine, light brown, medium brown, olive, red, blue, grey, illusion, gingiva; clear, clear h.c.; CREATION VC Translucent h.v., Translucent I.v.) | |
| 2. Nachfüllpackungen | |
| 15 Make up Neo stain A, stain B, stain C, white, eggshell, lemon yellow, mandarine, light brown, medium brown, olive, red, blue, grey, illusion, gingiva | 2,6 ml |
| 2 CREATION VC Translucent h.v., Translucent I.v. | 2,4 ml |
| 2 Make up Neo clear, clear h.c. | 5 ml |

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PLEASE READ CAREFULLY BEFORE USE. KEEP FOR REFERRAL PURPOSES.

EN

CREATION VC Make up Neo

Light-curing characterisation material.

For use only by a dental technician in the recommended indications.

INDICATIONS

- Characterization of direct and indirect composite restorations, acrylic denture bases and artificial acrylic teeth.
- For obtaining surface smoothness and wear resistance of restorations made of composite resin, acrylic denture bases and artificial acrylic teeth.

Note: not suitable for characterization of natural teeth.

KONTRAINDICATIONS

Avoid use of this product in patients with known allergies to acrylate / methacrylate monomers or acrylate / methacrylate polymers.

USE

Surface characterization

- After contouring the resin surface, finish with a carbide bur or coarse silicone point. The oxygen inhibition layer must be removed as it inhibits the polymerization of CREATION VC Make up Neo.

Laboratory use: sandblast with 25-50 µm aluminium oxide (0.15 MPa/1.5 bar), clean and dry.

Note:

- If surface still contains oil residue, the surface can be cleaned with ethanol (alcohol.) The use of ethanol (alcohol) on acrylic resin may cause microscopic cracks. A mild detergent can be used to clean the acrylic resin.
- When applying CREATION VC Make up Neo on hybrid ceramic blocks treat the surface using CREATION VC Ceramic Primer.
- Shake the bottle of CREATION VC Make up Neo well. Dispense a few drops into a dispensing dish. Dilute by adding CREATION VC Make up Neo clear to the desired shade. Apply a thin layer to the resin surface using the brush included in the delivery. Do not air blow.
- Light cure layer with a suitable light curing device (see chart 1).

Laboratory use: set a 2 cm platform in the curing device and place the prosthesis on it with the coated surface facing the light. Reverse the prosthesis and repeat the light curing to fully cure the undercut or area in shadow.

Internal characterization

1. Build up the dentine layer with CREATION VC Veneering Composite and light cure. Observe the instructions of use of CREATION VC Veneering Composite.

Note: do not remove the oxygen inhibition layer. Apply CREATION VC Make up Neo on the oxygen inhibition layer. If the oxygen inhibition layer was removed, treat the surface with CREATION VC Ceramic Primer (dry with compressed air) and CREATION VC Composite Primer (light cure) before applying CREATION VC Make up Neo.

- Shake the bottle of CREATION VC Make up Neo well. Dispense a few drops into a dispensing dish. Apply a thin layer of CREATION VC Make up Neo of choice to the dentine surface with an instrument or brush (e.g. the convex and concave of dentine mamelon or incisal area). Mix multiple colours to obtain the desired shade if needed. Do not dry with air.
- Light cure layer with a suitable light curing device (see chart 1).
- Continue building up with CREATION VC Enamel, Effect Enamel, Transpa on the oxygen inhibition layer. Observe the instructions of use of CREATION VC Veneering Composites.
- Finish and polish using standard technique or apply CREATION VC Make up Neo, clear or clear h.c. on the entire surface.

Note: for diluting or tinting, CREATION VC Translucent h.v. and Translucent I.v. are especially suitable for Make up Neo colours. These two materials are only suitable for internal characterisations.

Repair or restoration (e.g. when the gloss is lost)

- Roughen the surface to be treated using a carbide bur and paper.
 - Clean the surface with oil free air.
- Note:** if surface still contains oil residue, the surface can be cleaned with ethanol (alcohol.) The use of ethanol (alcohol) on acrylic resin may cause microscopic cracks.
- Apply CREATION VC Make up Neo, clear or clear h.c. shade to the surface.Do not dry with air.
 - Light cure with a suitable light curing device (see chart 1).

P LEASE READ CAREFULLY BEFORE USE. KEEP FOR REFERRAL PURPOSES.

EN

Chart 1: light curing time

| | | |
|---|-----------------------------------|---------|
| GC LABOLIGHT DUO* | | 90 sec. |
| GC LABOLIGHT LV-III, II* | | 5 min. |
| GC STEPLIGHT SL-1* | distance from light source < 3 cm | 10 sec. |
| | distance from light source < 3 cm | 20 sec. |
| Halogen light | | 40 sec. |
| Plasma arc | | 8 sec. |
| LED light (with light source having wavelength lower than 430 nm) | | 40 sec. |

Chart 2: depth of cure

| | |
|---|---------|
| MUN-4 mandarine, MUN-3 lemon yellow, MUN-8 medium brown, MUN-G gingiva | 0,07 mm |
| MUN-A stain A, MUN-B stain B, MUN-C stain C, MUN-1 white, MUN-2 eggshell, MUN-7 light brown, MUN-10 red, MUN-9 olive, MUN-ILilusion, MUN-12 grey, MUN-11 blue | 0,15 mm |
| MUN-CL clear, MUN-CL clear h.c. | 0,50 mm |

STORAGE

Recommended for optimal performance: store at room temperature of (4-25 °C).

SIZE OF PACKAGING

1. Set Packages

CREATION VC Make up Neo Kit (Make up Neo stain A, stain B, stain C, white, eggshell, lemon yellow, mandarine, light brown, medium brown, olive, red, blue, grey, illusion, gingiva; clear, clear h.c.; CREATION VC Translucent h.v., Translucent I.v.)

2. Refill Packages

- | | |
|--|--------|
| 15 Make up Neo stain A, stain B, stain C, white, eggshell, lemon yellow, mandarine, light brown, medium brown, olive, red, blue, grey, illusion, gingiva | 2,6 ml |
| 2 CREATION VC Translucent h.v., Translucent I.v. | 2,4 ml |
| 2 Make up Neo clear, clear h.c. | 5 ml |

CAUTION

- Always wear personal protective equipment (PPE) such as gloves, face masks and safety eyewear.
- In particular avoid contact of uncured material with skin.
- Avoid getting material on clothing.
- CREATION VC Make up Neo is highly flammable and volatile. Do not use near open flame. Keep away from sources of ignition.
- CREATION VC Make up Neo is volatile and may cause respiratory irritation. Use only in well ventilated rooms. Replace cap immediately after use.
- Do not use CREATION VC Make up Neo on high pressure bearing areas.
- In case of contact with eyes, immediately flush with water and seek medical attention.
- In case of contact with skin, immediately remove with a sponge or cotton pellet. Flush with water.
- Remove the oxygen inhibition layer of the composite surface to be coated before applying CREATION VC Make up Neo as it may hinder the proper curing of CREATION VC Make up Neo.
- Dispense a few drops of CREATION VC Make up Neo into a disposable dispensing dish immediately before applying as the material contains a volatile solvent.
- Store locked up.
- Dispose all waste according to local regulations.
- In rare cases the product may cause sensitivity to some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

CLEANING AND DISINFECTING MULTI-USE DELIVERY SYSTEMS

To avoid cross-contamination between patients, disinfection with suitable agents is required. Immediately after use inspect product for deterioration. Discard device if damaged.
D . Thoroughly clean material and protect against humidity to prevent cross-contamination. Disinfect with a (mid-level registered) healthcare-grade infection control product according to regional/national guidelines.

* Observe the manufacturer’s instructions for use and make sure that these devices have the same light output as the GC LABOLIGHT LV-II, III, GC LABOLIGHT DUO and GC STEPLIGHT SL-I.

Some products referenced in this IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: https://www.creation-willigeller.com/en/media-center/manuals-brochures/

Undesired effects - Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental. In this way you will contribute to improve the safety of this product.

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1. Assortiment

CREATION VC Make up Neo Kit (Make up Neo stain A, stain B, stain C, white, eggshell, lemon yellow, mandarine, light brown, medium brown, olive, red, blue, grey, illusion, gingiva; clear, clear h.c.; CREATION VC Translucent h.v., Translucent l.v.)

2. Réassorts

15 Make up Neo 2,6 ml
stain A, stain B, stain C, white, eggshell, lemon yellow, mandarine, light brown, medium brown, olive, red, blue, grey, illusion, gingiva
2 CREATION VC Translucent h.v., Translucent l.v. 2,4 ml
2 Make up Neo clear, clear h.c. 5 ml

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© Creation Willi Geller International GmbH 02/2021 - Rev. 5 | 829915
Distributed by: GC AUSTRALASIA DENTAL PTY. LTD., 1753 Botany Rd, Banksmeadow, NSW 2019, Australia, Tel. +61 2 9301 8200
Distributor: Creation Willi Geller International GmbH, Koblacherstraße 3, 6812 Meiningen, Austria