

WILLI GELLER
Creation
PRODUCT CATALOGUE



WILLI GELLER
Creation
MAGIC COLOUR
Highlight
white
HL-1


WILLI GELLER
Creation
MAGIC COLOUR
Highlight
eggshell
HL-2

WILLI GELLER
Creation
MAGIC COLOUR
Make-In
honey yellow
MI-61

WILLI GELLER
Creation
MAGIC COLOUR
DS-A

WILLI GELLER
Creation
MAGIC COLOUR
HL-1 white

2024



“I have no secret – my success is a result of all my experiences.”

Simply Brilliant – then as now

Inspired by the desire for a unique and vital copy of nature’s own creation, the internationally renowned dental technician Willi Geller created a metal and ceramics range in 1988 which made quite a splash: **perfect, almost natural restorations** became possible all at once – masterly in their physical properties and reliable and reproducible in terms of shading. With Creation, each tooth was suddenly able to retain its own character – in shade, form and function

The new ceramic materials from Switzerland quickly set the trend for dental technology worldwide. Geller’s **“master ceramics”** and his distinctive layering technique set benchmarks for a top-class dental aesthetic. Through the opportunity for ideal natural replicas, and through the seemingly **inexhaustible variety** of the metal and ceramics range, the dental technician gained an entirely new profile – the profile of an artist.

Experience as a source of success

Founded as a “one-man company” 30 years ago, the dental company has continued to grow and develop steadily: new manufacturing technologies, additional fusions and effect materials, as well as **expanded indications and, more recently, a light-curing veneer composite with “Creation Brilliance”**. Today, Creation is one of the world’s most important manufacturers of ambitious ceramics systems, and resin and composite teeth with a strong profile. The innovative ceramic and composite materials and denture teeth are successfully processed in over 70 countries **across the globe**. And the users are not the only ones who are impressed; dentists and patients are also convinced by our “Simply Brilliant“ aesthetics concept.

Willi Geller’s dental creations represent the original – and are just as brilliant today.

WILLI GELLER
Creation

Creation Highlights

VENEERING CERAMIC Page

Metal

Creation CC – Classic Metal Ceramic 05

Zirconium

Creation ZI-CT – Zirconium Dioxide Ceramic 10

Creation ZI-F – Low-Fusing Zirconium Dioxide Ceramic 14

Lithium Disilicate

Creation LS/LS Press – Lithium Disilicate
Veneering and Pressable Ceramic 18

OT+ AND GELLER-MIX Page

Opal Transpa OT+ and Geller-Mix 24

COLOURS Page

Magic Colour 26

CreaColor

Make up Neo 30

In Nova Neo 30

LIQUIDS Page

Modelling Liquid – High & Low Fusing 32

Opaque Liquid – High Fusing 32

Shoulder Liquid 32

Universal & Glaze Liquid 32

Moisturing Fluid 33

Magic Colour Liquid GL & Special Liquid ML 33

Make up Neo Fluid & Reconditioner 33

In Nova Neo Fluid 33

SHADE GUIDES Page

Creation Shade Guides 34

CERAMIC ACCESSORIES Page

CreaPen-Separating Pen 36

Liquid Firing Pad 36

Creation Brush Kit 36

Aqualine master 37

Aqualine smart 37


Aqualine mini 37

PRODUCT LIST Page

Product List 40



REM-picture,
etched, ceramic surface:



Leucite-reinforced metal ceramic CC

CREATION CC – technically refined – inspired by nature

Creation CC is a high-fusing metal ceramic with unique optical and physical properties that has developed constantly in line with growing aesthetic demands – **recognized worldwide for 30 years!**

The perfectly coordinated ceramic materials have excellent homogeneity and thus guarantee high flexural strength. The result: a densely sintered structure for **pure and non-porous layering** with the utmost reliability.

The pure potash feldspar with micro-pure leucite crystal structure in Creation CC creates iridescent light refraction and natural brilliance and enhances the strength. Exceptional colour nuances can additionally be created with the unique effect materials.

The ceramic bonder Crea Alloy Bond also acts as a CTE buffer between metal and ceramic – for improved adhesive properties to non-precious metal and greater colour reliability.

With the modular classic brand Creation CC, metal-ceramic restorations can be produced with impressive aesthetics and dimensional stability – for **superb dental artistry!**

Creation CC – natural brilliance and efficiency.

THE BRILLIANCE OF CREATION CC:

- Unique, lifelike aesthetics thanks to highly pure potash feldspar
- Resilience due to consistently high flexural strength
- Natural colour effect and light dynamics due to leucite crystals
- Ease of handling thanks to consistent colour and layering system
- Processing reliability since 1988

6 Creation CC – Classic Metal Ceramic



CC Creapast Kit | 000300
Complete paste opaques with modifier

- 16 Opaque O-A1-D4 3g
- 5 Opaque Dentine Modifier . 3g
- 1 Creapast Brush
- 1 CC Manual



CC Opaque Kit | 000800
Complete powder opaques with modifier

- 16 Opaque O-A1-D4 20g
- 8 Opaque Modifier
- OM-11-17, G 15g
- 1 Foundation Opaque
- FOP 20g
- 2 Opaque Liquid OL-H . . 50ml
- 1 Universal & Glaze Liquid 25ml
- 1 Shade Guide
- 1 CC Color Chart
- 1 CC Manual



CC Kit 1 | 000804
Basic porcelains kit (opaque dentines, dentines, incisals)

- 16 Opaque Dentine
- OD-A1-D4 15g
- 16 Dentine D-A1-D4 20g
- 4 Enamel E-57-60 20g
- 2 Transpa NT, OT 20g
- 1 Clear CL-O 20g
- 1 Glaze GL 10g
- 2 Modelling Liquid ML-H . 50ml
- 1 Universal & Glaze Liquid 25ml
- 2 Shade Guides
- 1 CC Color Chart
- 1 CC Manual



CC Kit 2 | 001508
Effect porcelains kit

- 5 Opaque Dentine
- OD-32, 37, 41, 43, 44 . . 15g
- 1 Dentine D-CW 20g
- 2 Approximal Dentine
- AD-1, 2 15g
- 1 Clear UC 20g
- 5 Transpa TI-1-5 20g
- 8 Effect Enamel
- SI-01-06, SO-10, 11 . . . 20g
- 4 Pearl Enamel PS-0-3 . . . 20g
- 6 Neck Transpa HT-51-56 15g
- 6 Make In MI-61-66 15g
- 1 Correction Powder KM . . 20g
- 2 Modelling Liquid ML-H 50ml
- 1 Universal & Glaze Liquid 25ml
- 2 Shade Guides
- 1 CC Color Chart
- 1 CC Manual



CC Bleach Shades Kit | 000250

Complete bleach shades colours

- | | |
|---|-----------------------------------|
| 1 Bleach Creapast Opaque
O-AB 3g | 1 Shoulder Liquid 25 ml |
| 3 Bleach Dentine
BD-A, B, BO 20g | 1 Modelling Liquid ML-H . 25 ml |
| 1 Bleach Enamel S-AB . . . 20g | 1 Bleach Shade Guide |
| 1 Bleach Shoulder SP-AB . 15g | 1 CC Color Chart/Firing Chart |



CC Gingiva Kit | 001590

Complete gingiva colours

- | | |
|---------------------------------|---------------------------------|
| 1 Creapast Opaque O-G . . . 3g | 1 Opaque Liquid OL-H . . 25 ml |
| 1 Opaque Modifier OM-G . 15g | 1 Modelling Liquid ML-H . 25 ml |
| 8 Gingiva G1-7, N 15g | 1 Shade Guide |



CC Shoulder Kit | 001503

Complete shoulder porcelains

- | |
|---|
| 10 Shoulder Powder
SP-21-29, G 15g |
| 1 Shoulder Liquid 25 ml |
| 1 CC Color Chart/Firing Chart |



CC Transition Dentine Kit | 000801

Basic porcelains kit (transition dentines, incisals)

- | | |
|---|----------------------------------|
| 16 Opaque Dentine
OD-A1-D4 15g | 2 Modelling Liquid ML-H . 50 ml |
| 16 Transition Dentine
TD-A1-D4 20g | 1 Universal & Glaze Liquid 25 ml |
| 4 Enamel E-57-60 20g | 1 Shade Guide |
| 1 Clear CL-O 20g | 1 CC Color Chart |
| 2 Transpa NT, OT 20g | 1 CC Manual |
| 1 Glaze GL 10g | |

8 Creation CC – Classic Metal Ceramic



CC Starter Kit | 000915

Introduction kit with tooth shades A2, A3, A3.5, B2, B3

- | | |
|---------------------------------|-----------------------------------|
| 2 Crea Alloy Bond 4g | 5 Transpa NT, OT, |
| 5 Creapast Opaque. 3g | TI-1, 4, 5 20g |
| 5 Opaque Dentine 15g | 2 Effect Enamel SI-02, 04 . 20g |
| 3 Opaque Dentine | 3 Neck Transpa HT-52-54 15g |
| OD-32, 37, 43 15g | 1 Glaze GL 10g |
| 5 Dentine 20g | 2 Make In MI-63, 65 15g |
| 2 Enamel E-58, 59 20g | 2 Modelling Liquid ML-H . 25ml |
| 1 Clear CL-O 20g | 2 Universal & Glaze Liquid 25ml |
| | 1 CC Manual |



CC Introduction Kit | 000920

Introduction kit with tooth shades A1, A2, A3, A3.5, B1, B2, C1, C2, C3, D2, D3

- | | |
|------------------------------------|-----------------------------------|
| 11 Creapast Opaque. 3g | 3 Effect Enamel |
| 3 Opaque Dentine | SI-02, 04, SO-10 20g |
| OD-32, 37, 43 15g | 1 Pearl Enamel PS-1 20g |
| 11 Dentine 20g | 3 Neck Transpa HT-52-54 15g |
| 4 Enamel E-57-60 20g | 1 Glaze GL 10g |
| 1 Clear CL-O 20g | 1 Modelling Liquid ML-H . 50ml |
| 4 Transpa NT, TI-1, 2, 4 . . . 20g | 1 Universal & Glaze Liquid 25ml |
| | 1 CC Color Chart & CC Manual |

Crea Alloy Bond | 000080

Crea Alloy Bond is a bonding agent which was specially developed for neutralising leaking metallic oxides. At the same time it acts as an efficient CTE buffer between porcelain and metal alloy.



- Safety owing to high flexural strength (70 MPa)
- No marginal discolorations – staining metallic oxides are neutralised
- Greater colour reliability because of the golden yellow undertone
- Compatible with all materials from Creation CC and LF

WOP Opaque

The coating of opaque can be carried out at only 950 °C instead of 980 °C with the WOP paste. This procedure is useful with certain high-gold alloys and galvanised frameworks. After the WOP coating the second coating of opaque is applied.



- Firing temperature 950 °C

- | |
|--|
| CC Creapast WOP sand 3g 000322 |
| CC Creapast WOP sahara 3g 000323 |
| CC Creapast WOP gold 3g 000324 |


A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
OPAQUE DENTINE															
A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
DENTINE															
TD-A1	TD-A2	TD-A3	TD-3,5	TD-A4	TD-B1	TD-B2	TD-B3	TD-B4	TD-C1	TD-C2	TD-C3	TD-C4	TD-D2	TD-D3	TD-D4
TRANSITION DENTINE															
S-57	S-58	S-59	S-60		OT+	OT opal	Geller-Mix	NT neutral	TI-01 blue	TI-02 white	TI-03 pink	TI-04 yellow	TI-05 grey		
ENAMEL				TRANSPA											
SI-01 light yellow	SI-02 medium yellow	SI-03 heavy yellow	SI-04 light orange	SI-05 medium orange	SI-06 heavy orange	SO-10 blue	SO-11 orange		PS-0 white	PS-1 blue/w.	PS-2 flamingo	PS-3 orange		CL-O clear	UC window
EFFECT ENAMEL						PEARL ENAMEL						CLEAR			
HT-51 iridescent	HT-52 khaki	HT-53 sand	HT-54 honey	HT-55 borde.	HT-56 ocher	SP-21 neutral	SP-22 flamingo	SP-23 sand	SP-24 gold	SP-25 red brown	SP-26 yellow gold	SP-27 ivory	SP-28 opaque yellow	SP-29 honey	SP-G gingival
NECK TRANSPA HIGH FLUORESCENT						SHOULDER POWDERS HIGH FLUORESCENT									
MI-61 ivory	MI-62 lemon	MI-63 honey yellow	MI-64 flamingo	MI-65 gold	MI-66 olive		AD-1 light yellow	AD-2 orange							
MAKE IN HIGH FLUORESCENT INTERNAL POWDER						APPROXIMAL DENTINE									
G1 purple	G2 dark pink	G3 light pink	G4 flamingo	G5 rose	G6 dark pink opaque	G7 violet	G-N neutral		OD-32 havanna	OD-37 curry	OD-41 orange	OD-43 ivory	OD-44 cuba		
GINGIVAL						OPAQUE DENTINE									
O-AB BLEACH OPAQUE	BD-A BLEACH DENTINE	BD-B	BD-BO		TD-BA TRANSITION DENTINE	S-AB BLEACH ENAMEL		SP-AB BLEACH SHOULDER		KM CORRECTION POWDER				GL GLAZE	

Shades A1–D4 are referring to the Vita® Shade Guide. Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.


10 Creation ZI-CT – Zirconium Dioxide Ceramic




REM-pictures
etched, ceramic surfaces:



Glass ceramic



Leucite-reinforced metal ceramic CC



Leucite-reinforced
zirconium dioxide ceramics ZI-CT

CREATION ZI-CT – the fine art of ceramics

Creation ZI-CT is an innovative feldspar ceramic containing leucite crystals which is matched to the CTE of all high-strength zirconium dioxide frameworks and structures – whether white, coloured or translucent.

The high proportion of pure potash feldspars creates light scatter that is comparable to that of natural teeth and a **brilliant colour effect**, while the micro-fine leucite crystals improve mechanical strength and thereby ensure a **durable ceramic structure**.

Fascinatingly fluorescent: the special bonding agent Frame Shade NT guarantees strong adhesion to the zirconium dioxide framework and also increases its fluorescence – for **powerful bond strength** without tensile stresses.

The modular Creation ZI-CT allows restorative work on zirconium dioxide with impressive colour and dimensional stability and with **brilliant light reflection and dynamics**.

High-fusing, high aesthetics, high stability: Creation ZI-CT.

“Unobtrusively beautiful – a true Creation!”

Andreas Nolte, MDT

THE BRILLIANCE OF CREATION ZI-CT

- Natural colour effect and light dynamics due to leucite crystals
- Excellent bonding strength thanks to high firing temperature of 910 °C
- Extremely low shrinkage for fewer correction firings and more efficiency
- Reliability due to homogeneous crystalline structure and high resilience



ZI-CT Kit 1 | 620501

Basic porcelains kit (dentines, incisors)

- | | |
|--|---------------------------------|
| 1 Frame Shade NT 4g | 1 Correction Powder KM . 20g |
| 5 Opaque Dentine
OD-32, 37, 41, 43, 44. . . 15g | 1 Glaze GL 10g |
| 16 Dentine D-A1-D4. 20g | 1 Universal & Glaze Liquid 25ml |
| 4 Enamel E-57-60 20g | 2 Modelling Liquid ML-H . 50ml |
| 1 Clear CL-O 20g | 1 ZI-CT Color Chart & Manual |
| 2 Transpa NT, OT 20g | |
| 2 Approximal Dentine
AD-1, 2. 15g | |



ZI-CT Kit 2 | 620502

Effect porcelains kit

- | | |
|---|---------------------------------|
| 1 Clear UC 20g | 1 Universal & Glaze Liquid 25ml |
| 5 Transpa TI-1-5. 20g | 2 Modelling Liquid ML-H . 50ml |
| 8 Effect Enamel
SI-01-06, SO-10, 11. . . 20g | 1 ZI-CT Color Chart & Manual |
| 4 Pearl Enamel PS-0-3. . . 20g | |
| 6 Neck Transpa
HT-51-56. 15g | |
| 6 Make In MI-61-66 15g | |



ZI-CT Starter Kit | 620503

Introduction kit with tooth shades A1, A2, A3, A3.5, B2, B3

- | | |
|---|----------------------------------|
| 1 Frame Shade NT 4g | 1 Pearl Enamel PS-0. 20g |
| 3 Opaque Dentine
OD-32, 37, 43 15g | 3 Neck Transpa HT-52-54 15g |
| 6 Dentine 20g | 3 Shoulder SP-21, 23, 28 . 15g |
| 2 Enamel E-58, 59 20g | 2 Make In MI-63, 65 15g |
| 1 Clear CL-O 20g | 1 Glaze GL 10g |
| 6 Transpa NT, OT
TI-1, 3-5 20g | 2 Modelling Liquid ML-H . 25ml |
| 2 Effect Enamel SI-02, 04 20g | 1 Shoulder Liquid 25ml |
| | 1 Universal & Glaze Liquid 25ml |
| | 1 ZI-CT Color Chart & Manual |



ZI-CT Introduction Kit | 620506

Introduction kit with tooth shades A1, A2, A3, A3.5, B1, B2, C1, C2, C3, D2, D3

- | | |
|---|---|
| 1 Frame Shade NT 4g | 1 Pearl Enamel PS-0. 20g |
| 3 Opaque Dentine
OD-32, 37, 43 15g | 3 Neck Transpa
HT-52-54. 15g |
| 11 Dentine 20g | 1 Glaze GL 10g |
| 1 Clear CL-O 20g | 2 Modelling Liquid ML-H . 50ml |
| 4 Enamel E-57-60 20g | 1 Universal & Glaze Liquid 25ml |
| 4 Transpa NT, TI-1, 3, 4 . . 20g | 1 ZI-CT Color Chart & Manual |
| 1 Effect Enamel SO-10 . . . 20g | |

12 Creation ZI-CT – Zirconium Dioxide Ceramic



ZI-CT Bleach Shades Kit | 620505
Complete bleach shades colours

- 3 Bleach Dentine
BD-A, B, BO 20g
- 1 Bleach Enamel S-AB . . . 20g
- 1 Shoulder SP-AB. 15g
- 1 Modelling Liquid ML-H 25ml
- 1 Shoulder Liquid 25ml
- 1 ZI-CT Color Chart



ZI-CT Gingiva Kit | 620504
Complete gingiva colours

- 8 Gingiva G1–7, N. 15g
- 1 Modelling Liquid ML-H 25ml
- 1 ZI-CT Color Chart
Bleach & Gingiva



ZI-CT Shoulder Kit | 620507
Complete shoulder porcelains

- 10 Shoulder
SP-21–29, G 15g
 - 1 Shoulder Liquid 25ml
 - 1 ZI-CT Color Chart/Firing Chart
-

A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
DENTINE															
E-57	E-58	E-59	E-60		OT+	OT opal	Geller-Mix	NT neutral	TI-01 blue	TI-02 white	TI-03 pink	TI-04 yellow	TI-05 grey		
ENAMEL				TRANSPA											
SI-01 light yellow	SI-02 medium yellow	SI-03 heavy yellow	SI-04 light orange	SI-05 medium orange	SI-06 heavy orange	SO-10 blue	SO-11 orange		PS-0 white	PS-1 blue/w.	PS-2 flamingo	PS-3 orange		CL-O clear	UC window
EFFECT ENAMEL						PEARL ENAMEL						CLEAR			
HT-51 iridescent	HT-52 khaki	HT-53 sand	HT-54 honey	HT-55 borde.	HT-56 ocher	SP-21 neutral	SP-22 flamingo	SP-23 sand	SP-24 gold	SP-25 red brown	SP-26 yellow gold	SP-27 ivory	SP-28 opaque yellow	SP-29 honey	
NECK TRANSPA HIGH FLUORESCENT						SHOULDER POWDERS HIGH FLUORESCENT									
MI-61 ivory	MI-62 lemon	MI-63 honey yellow	MI-64 flamingo	MI-65 gold	MI-66 olive		AD-1 light yellow	AD-2 orange							
MAKE IN HIGH FLUORESCENT INTERNAL POWDER						APPROXIMAL DENTINE									
G1 purple	G2 dark pink	G3 light pink	G4 flamingo	G5 rose	G6 dark pink opaque	G7 violet	G-N neutral		OD-32 havanna	OD-37 curry	OD-41 orange	OD-43 ivory	OD-44 cuba		
GINGIVAL								OPAQUE DENTINE							
BD-A	BD-B	BD-BO		S-AB		SP-AB		KM							GL
BLEACH DENTINE			BLEACH EN.			BLEACH SHOULDER			CORRECTION POWDER					GLAZE	

Shades A1–D4 are referring to the Vita® Shade Guide. Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

14 Creation ZI-F – Low-Fusing Zirconium Dioxide Ceramic



CREATION ZI-F – naturally brilliant – brilliantly natural

Creation ZI-F is a tried and tested zirconium dioxide ceramic. It is perfectly suited to veneering any kind of zirconium dioxide framework.

At the same time, the modular materials are also an ideal choice as **add-on materials** for Creation ZI-CT because their firing temperature is about 100°C lower than the high-fusing zirconium dioxide veneering ceramic.

Creation ZI-F is an **extremely versatile** zirconium dioxide ceramic with high **natural light reflection**, combined with all the material-related advantages of zirconium dioxide.

With the unique effect materials, exceptional colour nuances can also be created rapidly and comfortably – for maximum impact with minimum effort.

Creation ZI-F – natural light dynamics!



ZI-F Kit 1 | 653010

Basic porcelains kit (dentines, incisals)

- | | |
|---|----------------------------------|
| 5 Frame Shade | 2 Modelling Liquid ML-L . 50 ml |
| FS-1-4, IL 10 g | 1 Universal & Glaze Liquid 25 ml |
| 16 Dentine D-A1-D4 20 g | 2 Shade Guides |
| 4 Enamel E-57-60 20 g | 1 ZI-F Color Chart |
| 2 Clear CL-O, UC 20 g | 1 ZI-F Manual |
| 6 Transpa NT, OT,
TI-1, 3-5 20 g | |
| 1 Glaze GL 10 g | |



ZI-F Kit 2 | 653015

Effect porcelains kit

- | | |
|----------------------------------|-----------------------------------|
| 5 Opaque Dentine | 2 Approximal Dentine |
| OD-32, 37, 41, 43, 44 . . 15 g | AD-1, 2. 15 g |
| 5 Effect Enamel | 1 Shoulder Liquid 25 ml |
| SI-02, 04, 06, SO-10, 11 20 g | 2 Modelling Liquid ML-L . 25 ml |
| 2 Pearl Enamel PS-0, 3 . . 20 g | 1 Universal & Glaze Liquid 25 ml |
| 3 Neck Transpa HT-52-54 15 g | 2 Shade Guides |
| 5 Shoulder SP-22-25, 28 . 15 g | 1 ZI-F Color Chart |
| 3 Make In MI-61, 63, 65 . . 15 g | 1 ZI-F Manual |
| | 1 ZI-F Color Chart/Firing Chart |



ZI-F Bleach Shades Kit | 652050

Complete bleach shades colours

- | | |
|---|-----------------------------------|
| 1 Bleach Frame Shade
FS-AB 10g | 1 Shoulder Liquid 25 ml |
| 3 Bleach Dentine
BD-A, B, BO 20g | 1 Modelling Liquid ML-L . 25 ml |
| 1 Bleach Enamel S-AB . . . 20g | 1 Shade Guide |
| 1 Bleach Shoulder SP-AB . 15g | 1 ZI-F Color Chart/Firing Chart |



ZI-F Gingiva Kit | 603030

Complete gingiva colours

- | |
|---------------------------------|
| 1 Frame Shade FS-G 10g |
| 8 Gingiva G1-7, N 15g |
| 1 Modelling Liquid ML-L . 25 ml |
| 1 Shade Guide |



ZI-F Starter Kit | 653025

Introduction kit with tooth shades A1, A2, A3, A3.5, B2, B3


- | | |
|---|-----------------------------------|
| 2 Frame Shade FS-1, 2 . . . 10g | 3 Neck Transpa HT-52-54 15g |
| 3 Opaque Dentine
OD-32, 37, 43 15g | 3 Shoulder SP-22, 23, 28 . 15g |
| 6 Dentine 20g | 2 Make In MI-63, 65 15g |
| 2 Enamel E-58, 59 20g | 1 Glaze GL 10g |
| 1 Clear CL-O 20g | 1 Shoulder Liquid 25 ml |
| 6 Transpa NT, OT, TI-1, 3-5 20g | 2 Modelling Liquid ML-L . 25 ml |
| 2 Effect Enamel SI-02, 04 . 20g | 1 Universal & Glaze Liquid 25 ml |
| 1 Pearl Enamel PS-0 20g | 1 ZI-F Color Chart |
| | 1 ZI-F Manual |

ZI-F Single Products



- | | |
|-----------------------|----------------------|
| Frameshade FS-NT | 601087 4g |
| Shoulder Powder SP-21 | 602121 15g |
| Transpa TI-2 | 652077 20g |
| Neck Transpa HT-51 | 652151 15g |
| Neck Transpa HT-56 | 652156 15g |
| Correction Powder KM | 652090 20g |

16 Creation ZI-F – Color Chart



LS-Press
MT-B00



LS-Press
MO-1

LS-Press
WT-B00

LS-Press
LT-A





CREATION LS/LS Press – Press, veneer and individualise with lithium disilicate!

Efficiency and natural aesthetics: Creation LS Press and Creation LS are a perfect ceramic symbiosis offering dental technicians reliability and stability but also brilliant colours and versatility on and with lithium disilicate, as standard is not enough for Creation!

Creation Willi Geller completes the product range with a new pressable and veneering ceramic allowing the dental technician outstanding aesthetics, whether used single or in combination. With the Creation LS Press blanks and the Creation LS veneering ceramic perfectly matched to the lithium disilicate frame highly aesthetic and stable all-ceramic restorations can be manufactured – efficiently pressed and then individually veneered resp. characterised with the proven CreaColor stains and modifiers.

Reliably combined for extended possibilities: the two system components Creation LS Press and Creation LS follow master ceramist Willi Geller's "Simply Brilliant" philosophy proven for 30 years. They allow the dental technician manufacturing naturally beautiful anterior and posterior restorations with resp. on lithium disilicate assuring system safety.

Creation LS/LS Press – the sophisticated combination

THE BRILLIANCE OF THE CREATION LS:

- Comprehensive colour assortment for individual veneering and characterisation
- Layer build-up according to the tried-and-tested colour coding concept of Creation Willi Geller
- Natural colour brilliance and light effects due to feldspar content
- Homogeneous crystal structure and outstanding surface density
- Maximum stability – even after several firings

THE BRILLIANCE OF THE CREATION LS PRESS:

- High strength due to homogeneous crystalline structure: > 500 MPa
- Minimal invasive preparation: 0.3 mm veneer, 0.8 mm anterior
- Creative freedom of manufacturing: full veneers, partial veneers
- The appropriate blank for every indication

20 Creation LS/LS Press – Veneering / Pressable Ceramic



LS Kit | 830083

Basic and effect porcelains with the tooth shades BD-A, BD-B, BD-B0, A1, A2, A3, A3.5, B1, B2, B3, C2, D2

- | | |
|--|--|
| 4 Opaque Dentine
OD-32, 37, 41, 43 15 g | 3 Effect Enamel
SI-02, 04, SO-10 20 g |
| 12 Dentine BD-A, BD-B, BD-B0,
A1–A3.5, B1–B3, C2, D2 20 g | 1 Pearl Enamel PS-0 20 g |
| 4 Enamel E-57–60 20 g | 4 Neck Transpa HT-51–54. 15 g |
| 1 Bleach Enamel S-AB . . . 20 g | 2 Make-in MI-61, 63 15 g |
| 1 Clear CL-0 20 g | 1 Glaze GL 10 g |
| 7 Transpa NT, OT,
TI-01–05 20 g | 2 Modelling Liquid ML-L . 50 ml |
| | 1 Universal & Glaze Liquid 25 ml |
| | 1 Creation LS/LS Press Manual |



LS Press Starter Kit | 830080

Introduction kit with the tooth shades BD-A, BD-B0, A1, A2, A3 and matching press blanks with different opacities MO, LT, MT, HT

- | | |
|--|--|
| 3 LS Press MO-0–2 5 x 3 g | 1 Clear CL-0 20 g |
| 1 LS Press LT-A 5 x 3 g | 1 Transpa OT 20 g |
| 1 LS Press MT-B00 5 x 3 g | 2 Effect Enamel
SI-02, SO-10 20 g |
| 2 LS Press HT-2, 3 5 x 3 g | 1 Pearl Enamel PS-0 20 g |
| 2 Opaque Dentine
OD-32, 41 15 g | 1 Make-in MI-61, 63 15 g |
| 5 Dentine BD-A, BD-B0,
D-A1–A3 20 g | 1 Glaze GL 10 g |
| 2 Enamel E-58, 59 20 g | 1 Modelling Liquid ML-L . 50 ml |
| 1 Bleach Enamel S-AB . . . 20 g | 1 Universal & Glaze Liquid 25 ml |
| | 1 Creation LS/LS Press Manual |



LS Press Introduction Kit A1, A2 | 830081

The easy introduction into the system with the tooth shades A1, A2 and the matching press blanks with different opacities MO-1, LT-A, HT-2

- | | |
|-----------------------------------|----------------------------------|
| 1 LS Press MO-1 5 x 3 g | 1 Clear CL-0 20 g |
| 1 LS Press LT-A 5 x 3 g | 1 Transpa OT 20 g |
| 1 LS Press HT-2 5 x 3 g | 1 Glaze GL 10 g |
| 1 Opaque Dentine OD-32 . 15 g | 1 Modelling Liquid ML-L . 50 ml |
| 2 Dentine D-A1, A2 20 g | 1 Universal & Glaze Liquid 25 ml |
| 1 Enamel E-58 20 g | 1 Creation LS/LS Press Manual |



LS Press Introduction Kit A3 | 830082

The easy introduction into the system with the tooth shade A3 and the matching press blanks with different opacities MO-2, LT-A, HT-3

- | | |
|-----------------------------------|----------------------------------|
| 1 LS Press MO-2 5 x 3 g | 1 Clear CL-0 20 g |
| 1 LS Press LT-A 5 x 3 g | 1 Transpa OT 20 g |
| 1 LS Press HT-3 5 x 3 g | 1 Glaze GL 10 g |
| 1 Opaque Dentine OD-32 . 15 g | 1 Modelling Liquid ML-L . 50 ml |
| 2 Dentine D-A3 20 g | 1 Universal & Glaze Liquid 25 ml |
| 1 Enamel E-59 20 g | 1 Creation LS/LS Press Manual |
| 1 Effect Enamel SO-10 . . . 20 g | |



Creation LS Press MO-Ingot (Medium Opacity)

MO press ingots are the first choice for full veneering. Due to their high opacity, they feature a high light reflection and are suitable for covering discoloured dies and titanium abutments.

Colours

MO-0, MO-1, MO-2 5 x 3g



Creation LS Press LT-Ingot (Low Translucency)

LT press ingots can be used for restorations to be fully or partially veneered with Creation LS. With thin layer thicknesses, Creation LT can also be used for monolithic restorations.

Colours

LT-B00, LT-B0, LT-A1, LT-A2, LT-A3, LT-B1,
LT-A, LT-B, LT-C, LT-D 5x3g



Creation LS Press MT-Ingot (Medium Translucency)

MT press ingots are the first choice for fully-anatomic restorations which are characterised by individual staining. The MT ingot can also be used for full or partial veneering.

Colours

MT-B00, MT-B0, MT-A1, MT-A2, MT-A3,
MT-B1, MT-B2, MT-D2 5x3g



Creation LS Press HT-Ingot (High Translucency)

HT press ingots are particularly suitable for minimal invasive restorations (e.g. inlays, onlays and veneers) where mainly the natural enamel is replaced. The shade individualisation is carried out in the staining technique.

Colours

HT-B00, HT-B0, HT-1, HT-2, HT-3, HT-4 5x3g


22 Creation LS/LS Press – Veneering/Pressable Ceramic

Processing techniques


Opacity	Indications					Processing techniques		
	Occlusal veneers	Veneers	Inlays	Onlays	Crown	Staining technique	Cut-back technique	Layering technique
MO Medium Opacity					✓			✓
LT Low Translucency					✓		✓	✓
MT Medium Translucency	✓	✓	✓	✓	✓	✓	✓	
HT High Translucency	✓	✓	✓	✓	✓	✓		




Layering technique



Cut-back technique



Staining technique



Colour reference table

Tooth colour Vita®-Shade	A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
MO Medium Opacity	MO-1		MO-2			MO-1		MO-2		MO-1		MO-2				
LT Low Translucency	LT-A LT-A1	LT-A LT-A2	LT-A LT-A2/-A3		LT-A LT-A3	LT-B0 LT-B1	LT-B	LT-B	LT-B LT-A3	LT-C LT-B1	LT-C LT-A1	LT-C	LT-C	LT-D	LT-D	LT-D
MT Medium Translucency	MT-A1	MT-A2	MT-A2/-A3		MT-A3	MT-B0 MT-B1	MT-B2	MT-B2	MT-A3	MT-B1	MT-A1	MT-D2	MT-D2	MT-D2	MT-D2	MT-A2
HT High Translucency	HT-2		HT-3		HT-4	HT-1	HT-3		HT-4	HT-3	HT-4		HT-3			

Materials which are not related to the Vita®shade guide:

Bleach	MO-0	LT-B00	LT-B0	MT-B00	MT-B0	HT-B00	HT-B0	
--------	------	--------	-------	--------	-------	--------	-------	--

A1	A2	A3	A3,5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
DENTINE															
BD-A	BD-B	BD-BO		OD-32 havanna	OD-37 curry	OD-41 orange	OD-43 ivory		E-57	E-58	E-59	E-60		S-AB	
BLEACH DENTINE				OPAQUE DENTINE					ENAMEL				BLEACH EN.		
OT+	OT opal	Geller- Mix	NT neutral	TI-01 blue	TI-02 white	TI-03 pink	TI-04 yellow	TI-05 grey		SI-02 medium yellow	SI-04 light orange	SO-10 blue		PS-0 white	
TRANSPA										EFFECT ENAMEL			PEARL ENAMEL		
HT-51 irides- cent	HT-52 khaki	HT-53 sand	HT-54 honey		CL-O clear		MI-61 ivory	MI-63 honey yellow		KM				GL	
NECK TRANSPA HIGH FLUORESCENT					CLEAR		MAKE IN HIGH FLUORESCENT INTERNAL POWDER			CORRECTION POWDER			GLAZE		

LS-PRESS COLOR CHART

MO-0	MO-1	MO-2							
MEDIUM OPACITY									
LT-B00	LT-B0	LT-A1	LT-A2	LT-A3	LT-B1	LT-A	LT-B	LT-C	LT-D
LOW TRANSLUCENCY									
MT-B00	MT-B0	MT-A1	MT-A2	MT-A3	MT-B1	MT-B2	MT-D2		
MEDIUM TRANSLUCENCY									
HT-B00	HT-B0	HT-1	HT-2	HT-3	HT-4				
HIGH TRANSLUCENCY									

Shades A1–D4 are referring to the Vita® Shade Guide. Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.



Opal Transpa OT+ and Geller-Mix – highest opalescence and transparency for more liveliness!

With the new Geller-Mix and the Opal Transpa OT+, Creation Willi Geller expands its ceramic range with two layering ceramics with opalescent properties. The opalescence caused by light scattering and reflection increases the light dynamics of dental restorations and makes dentures appear even more natural in the incisal area.

The Geller-Mix and OT+ expand the range of all layering ceramics.

Opal Transpa OT+ and Geller-Mix – expressive in the colour effect and brilliant in the result!

- **Opal Transpa OT+** is a layering ceramic for maximum opalescence raising the light dynamic effects of veneered restorations to a higher level. OT+ is almost colour neutral and intensifies the opalescent light effects of both light and dark tooth shades.
- The **Geller-Mix** has been mixed individually by professionals for many years and is now available as a finished product. The fired ceramic with its subtle opacity convinces with its harmonious balance between opalescence and brightness. This property creates vivid effects at light ridges and incisal edges especially with light tooth shades.


Opal Transpa OT+

The effect porcelain OT+ shows intense opal effects and is suitable for both light and dark tooth shades.

Scope of application:
mesial and distal light ridges and incisal edges.



CC/Transpa OT+	20g 000287
CC/Transpa OT+	50g 000587
ZI-CT/Transpa OT+	20g 620046
ZI-CT/Transpa OT+	50g 620246
ZI-F/Transpa OT+	20g 652087
ZI-F/Transpa OT+	50g 655087
LS/Transpa OT+	20g 830056



Geller-Mix

The Geller-Mix is particularly used for light tooth shades due to its subtle milky opacity.

Scope of application:
mesial and distal light ridges and incisal edges.



CC/Transpa Geller-Mix	20 g	000288
CC/Transpa Geller-Mix	50 g	000588
ZI-CT/Transpa Geller-Mix	20 g	620047
ZI-CT/Transpa Geller-Mix	50 g	620247
ZI-F/Transpa Geller-Mix	20 g	652088
ZI-F/Transpa Geller-Mix	50 g	655088
LS/Transpa Geller-Mix	20 g	830057



CREATION Magic Colour – a colour system designed for the handcraft

Digital production of monolithic restorations is efficient and the results are robust. Magic Colour was also developed to meet users' aesthetic demands. The compact colour system connects innovative manufacturing processes with every individual dental technician's craftsmanship expertise.



Uncompromising aesthetics

Magic Colour is used for accentuating and surface staining of monolithic zirconia or lithium disilicate restorations. As a powder-liquid combination, the colour system opens up numerous processing options and is very easy to use: Mixed with Magic Colour Liquid GL, the finely ground colour powders are extremely efficient as paste compounds, even in thin layers, and the results are foreseeable.

More brilliance, more natural depth

Whether used for surface varnishing, for characterising colour effects, or in combination with microlayering: Magic Colour is the starting point for multifaceted aesthetics. Microlayering with Creation ZI-F (with zirconia) or with Creation LS (with lithium disilicate) enhances aesthetics, depth effect and brilliance – especially with monolithic anterior restorations. Magic Colour makes every restoration bear the signature of the respective dental technician.

THE BRILLIANCE OF MAGIC COLOUR:

- Manifold processing options thanks to the powder-liquid system
- Surface characterisation of CAD, press and layered ceramics, especially of monolithic zirconia and lithium disilicate restorations
- Tinting of layering ceramics
- For surface staining and for use in combination with Microlayering
- Natural, optical effect from the fluorescence of the fine-grained powders

28 Magic Colour



Magic Colour Professional Kit | 001205

The kit contains the full range of stains as well as glazes and liquids – for both white and red aesthetics.

- | | |
|-----------------------------------|------------------------------------|
| 4 Dentine Shade DS-A-D . . . 3g | 1 Special Liquid ML 50ml |
| 15 Highlight HL-1-15 3g | 1 Brush |
| 2 Make In MI-61, 63 3g | |
| 3 Gingiva Shade GS-1-3. . . 3g | |
| 1 Glaze 3g | |
| 1 Glaze Fluo 3g | |
| 1 Liquid GL 50ml | |



Magic Colour Starter Kit ZI-F | 001206

The Kit contains the most common stains, glazes and liquids for white aesthetics as well as two Creation ZI-F layering materials for micro layering zirconia restorations.

- | | |
|---|------------------------------------|
| 2 Dentine Shade DS-A, B . . 3g | 1 Glaze 3g |
| 6 Highlight HL-1, 4, 6, 8, 9, 11 3g | 1 Glaze Fluo 3g |
| 1 Make In MI-63 3g | 1 Liquid GL 25ml |
| 1 ZI-F Transpa NT 20g | 1 Special Liquid ML 25ml |
| 1 ZI-F Transpa OT 20g | 1 Brush |



Magic Colour Starter Kit LS | 001207

The Kit contains the most common stains, glazes and liquids for white aesthetics as well as two Creation LS layering materials for micro layering lithium disilicate restorations.

- | | |
|---|------------------------------------|
| 2 Dentine Shade DS-A, B . . 3g | 1 Glaze 3g |
| 6 Highlight HL-1, 4, 6, 8, 9, 11 3g | 1 Glaze Fluo 3g |
| 1 Make In MI-63 3g | 1 Liquid GL 25ml |
| 1 LS Transpa NT 20g | 1 Special Liquid ML 25ml |
| 1 LS Transpa OT 20g | 1 Brush |



created by Stefan Picha, Oral Design

MAGIC COLOUR CHART

DENTINE SHADE

DS-A	DS-B	DS-C	DS-D
------	------	------	------

HIGHLIGHT

HL-1 white	HL-2 eggshell	HL-3 lemon yellow	HL-4 mandar- ine	HL-5 flamingo	HL-6 apricot	HL-7 light brown	HL-8 medium brown	HL-9 olive	HL-10 red	HL-11 dove blue	HL-12 grey	HL-13 deep blue	HL-14 fissure	HL-15 illusion
---------------	------------------	-------------------------	------------------------	------------------	-----------------	------------------------	-------------------------	---------------	--------------	-----------------------	---------------	-----------------------	------------------	-------------------

MAKE IN

MI-61 ivory	MI-63 honey yellow
----------------	--------------------------

GINGIVA SHADE

GS-1 rasp- berry	GS-2 salmon red	GS-3 red violet
------------------------	-----------------------	--------------------

GLAZE

Glaze Fluo	Glaze
---------------	-------

LIQUIDE

Liquid GL 25ml	Liquid GL 50ml
----------------------	----------------------

SPECIAL LIQUIDE

Special Liquid ML 25ml	Special Liquid ML 50ml
---------------------------------	---------------------------------



CREACOLOR – stains and modifiers for all ceramics

CreaColor is a complete **“two in one” colour system** which can be used to place individual colour nuances simply and efficiently **for all ceramic surfaces and materials** – whether by the internal or external staining technique.

The physical properties of the **high-quality stains and modifiers** are so perfectly matched to the ceramic surfaces and materials that they can be worked universally with all layering and pressed ceramics, pressed or fired dentine cores as well as zirconium dioxide frames – after all, **colour creativity should be unlimited!** Furthermore, the new stains are extremely fine-grained and fluorescent and they can readily be mixed with each other.

CreaColor is made up of two components:

- CreaColor Make up Neo
- CreaColor In Nova Neo

THE BRILLIANCE OF CREACOLOR:

- Utmost colour brilliance for internal and external colouring
- Versatile and efficient: one system for all ceramics
- Natural colour effect due to very high fluorescence values
- Easy to handle: can be mixed with each other and can be corrected



CreaColor Kit | 001200
Contains Make up Neo Kit & In Nova Neo Kit



Make up Neo Kit | 001100
Glazing and fluorescent stains

- 13 Fluorescent Stain Glazes
 - MUN-F1–F13 4g
- 4 Dentine Shades MUN-FA–FD
in the four Vita® colour groups for colour corrections
and for painting onto press ingots 4g
- 1 Fluorescent Glaze GL-F 4g
- 1 Make up Neo Fluid 25ml
- 1 Make up Neo Reconditioner Fluid 8ml
(for freshening up concentrated stains)
- 1 Brush No. IV



In Nova Neo Kit | 001150
All-purpose modifiers

- 10 Fluorescent Modifiers
 - INN-1–10 3g
- 2 Effect Shades IL, CRL (Illusion and Crackliner)
for individual accents 3g
- 1 In Nova Neo Fluid 25ml
- 1 Brush No. IV

CREACOLOR COLOR CHART

MAKE UP NEO

MUN-F1 white	MUN-F2 eggshell	MUN-F3 lemon yellow	MUN-F4 mandarine	MUN-F5 flamingo	MUN-F6 apricot	MUN-F7 light brown	MUN-F8 medium brown	MUN-F9 olive	MUN-F10 red
MUN-F11 blue	MUN-F12 grey	MUN-F13 deep blue	MUN-FA stain A	MUN-FB stain B	MUN-FC stain C	MUN-FD stain D	MUN-GL-F Glaze	UF	Recon- ditioner

IN NOVA NEO

INN-1 melone	INN-2 rose	INN-3 yellow	INN-4 red brown	INN-5 brown	INN-6 olive	INN-7 grey	INN-8 blue	INN-9 white	INN-10 mais	INN-CRL Crackli- ner	INN-IL Illusion	INN- FLUID

WILLI GELLER
Creation

32 Creation Liquids

Fewer bottles – more simplicity and efficiency: For more clarity and space on the laboratory table, Creation liquids can be obtained in high fusion or low fusion versions.



Modelling Liquid – High & Low Fusing

Standard modelling liquid for mixing Creation ceramics except for shoulder ceramics (SP).

- ML-H (High) for high fusion ceramics Creation CC and ZI-CT
- ML-L (Low) for low fusion ceramic Creation ZI-F

Modelling Liquid High Fusing	25 ml 000012
Modelling Liquid High Fusing	50 ml 000013
Modelling Liquid High Fusing	250 ml 000014
Modelling Liquid High Fusing	1.000 ml 000015
Modelling Liquid Low Fusing	25 ml 000016
Modelling Liquid Low Fusing	50 ml 000017
Modelling Liquid Low Fusing	250 ml 000018



Opaque Liquid – High Fusing

Mixing liquid for the high fusion ceramic Creation CC.

Opaque Liquid High Fusing	25 ml 000019
Opaque Liquid High Fusing	50 ml 000020
Opaque Liquid High Fusing	250 ml 000021

Shoulder Liquid

Universal modelling liquid for mixing all SP shoulder ceramics from Creation.

Shoulder Liquid	25 ml 000025
Shoulder Liquid	50 ml 000026



Universal & Glaze Liquid

Additional universal liquid with a wide range of uses.

- Addition of a very small amount into the ceramic mixed with ML enhances the smoothness
- The pre-drying time must be extended during firing
- Mixing all Creation glazing materials and CreaColor colours In Nova Neo and Make up Neo
- Remoistening of Creapast

Universal & Glaze Liquid	25 ml 000010
Universal & Glaze Liquid	50 ml 000011



Moisturing Fluid

For extensive restorations and complex corrections, as a moisturiser for the applied ceramics.

- No more drying out, the ceramics remain consistently modellable
- Efficient and simple: one fluid for all ceramics
- Longer working time for mixed materials – especially for sizeable work
- Targeted application due to controllable moisture content

Moisturing Fluid 25 ml | 000090

Magic Colour Liquid GL

Mixing fluid for stain and glaze powder especially developed for the stain technique.

Liquid GL 25 ml | 001241

Liquid GL 50 ml | 001242



Magic Colour Special Liquid ML

Has been developed exclusively for microlayering with Creation ZI-F and Creation LS. It produces a homogeneous with dense ceramic consistency, which enables the application of material in very thin layers. Due to their low shrinkage and high density, the veneers exhibit a self-glaze effect after firing.

Special Liquid ML 25 ml | 001243

Special Liquid ML 50 ml | 001244

Make up Neo Fluid

Mixing liquid for Make up Neo.

Make up Neo Fluid 25 ml | 001125

Make up Neo Reconditioner

For freshening up concentrated stains.

Make up Neo Reconditioner 8 ml | 001126



In Nova Neo Fluid

Mixing liquid for In Nova Neo.

In Nova Neo Fluid 25 ml | 001163

KIT-BASED SHADE GUIDES

Art. no.	Article name
000064	CC/Shade Guide Kit 1 804
000073	CC/Shade Guide Kit 2 1508
000076	ZI-F/Shade Guide Kit 1 (28 Shades)
000077	ZI-F/Shade Guide Kit 2 (26 Shades)
000087	CC/Shade Guide Transition Dentine
000092	ZI-CT/Shade Guide Kit 1
000093	ZI-CT/Shade Guide Kit 2



SHADE GUIDES

Shades A1–D4 are referring to the Vita® Shade Guide.

Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.



Shade Guide Dentine A-D

Art. no. 000040




Shade Guide Opaque Dentine A-D

Art. no. 000041




Shade Guide Transpa

Art. no. 000042





Shade Guide Opaque Dentine + AD

Art. no. 000043





Shade Guide
Enamel

Art. no. 000044




Shade Guide
Transpa Modifier

Art. no. 000045





Shade Guide
Make In

Art. no. 000046





Shade Guide
Shoulder

Art. no. 000047




Shade Guide
Gingiva

Art. no. 000048




Shade Guide
Bleach

Art. no. 000049



Shade Guide
Wheel Empty

Art. no. 000091



Shade Guide
Single Shade Tab



Art. no. 000058 + Color description

36 Ceramic Accessories



CreaPen Separating Pen

The double-ended ceramic insulation pen with one pointed and one broad nib for surfaces or grooves respectively can be simply used on the sealed die.

CreaPen Separating Pen | 000059

Liquid Firing Pad

The liquid firing pad for all Creation systems precisely fixes the object on the fireproof tray.



- Stable hold on the firing tray
- Easy to handle because of proven dispenser design
- Safe because of firing stability and optimum heat/warmth build-up onto the object
- Economical because of precise measuring
- Can be used for all Creation systems
- Tough and easy to remove

Liquid Firing Pad | 000081

Creation Brush Kit

The six different sizes of the first-class Kolinsky red sable brush with a heavy wooden holder perfectly meet all the requirements of layering techniques.

Creation-Brush-Kit | 000136

1 Brush No. I | 000131

2 Brush No. II | 000132

3 Brush No. III | 000133

4 Brush No. IV | 000134

5 Brush No. V | 000135

6 Brush Creapast No. VI | 000321





aqua line



Aqualine Mixing Trays

With Aqualine, materials will always have the right viscosity. Drying out is prevented by the cellulose strips continuously absorbing liquid from the reservoir.

Aqualine master

Dimension: 335 x 225 x 25 mm

Components: Tray with glass plate (light grey) and hinged cover, self-adhesive write-on labels, cellulose strips

Aqualine master anthracite | AMA-001
Aqualine master white | AMA-002

Aqualine smart

Dimension: 260 x 170 x 32 mm

Components: Tray with glass plate (light grey) and hinged cover, self-adhesive write-on labels, cellulose strips

Aqualine smart anthracite | AS-001
Aqualine smart white | AS-002

Aqualine mini

Dimension: 133 x 71 x 23 mm

Components: Tray with glass plate (light grey) and hinged cover, cellulose strips

Aqualine mini anthracite | AM-001

ACCESSORIES

Aqualine master Glass Plate, printed light grey. 1 piece | AP-010

Aqualine master Cellulose Strips 40 pieces AP-013

Aqualine master Labels, self-adhesive for stain 5 pieces | AMA-018
for ceramic 30 pieces | AMA-019

Aqualine smart Tray Set Zircon stain plate 3 pieces ST-003

Aqualine master Tray Set Zircon stain plate 4 pieces ST-004

Aqualine smart Glass Plate, printed light grey. 1 piece | AS-010

Aqualine smart Cellulose Strips 30 pieces AS-012

Aqualine mini Tray Box with stain tray, with 15 cavities 1 piece AM-002

Aqualine mini Cellulose Strips 40 pieces AM-010

Aqualine Pool Set Dimension: 105 x 65 x 48 mm, rinsing cup 250 ml, sponge, sponge holder transparent. AS-014

Aqualine Pool Set Sponge 5 pieces AS-015

Aqualine Glass Care 2 x 30 ml AC-001

Aqualine Two in One Build up Liquid 500 ml AC-002



Creapearl 2

WILLI GELLER
Creation
SIMPLY BRILLIANT

WILLIGELLER
Creation
SIMPLY BRILLIANT

CONTENT

CREATION CERAMICS AND ACCESSORIES

All kits at a glance	41
CC - Classic Metal Ceramic	42
ZI-CT - Zirconium Dioxide Ceramic	51
ZI-F - Low-Fusing Zirconium Dioxide Ceramic	55
LS/LS Press - Lithium Disilicate Veneering and Pressable Ceramic	59
Magic Colour – Colour System	61
CreaColor – Stains and Modifier	62
Creation Liquids	63
Ceramic Accessories	63
Creation Shade Guides	64
Aqualine	65

Article no.	Article name	Content
	CERAMIC	
000300	CC/Kit Creapast	1 Kit
000800	CC/Powder Opaque Kit (inkl. Shade Guide 000068)	1 Kit
000804	CC/Kit 1 804 (inkl. Shade Guide 000064)	1 Kit
001508	CC/Kit 2 1508 (inkl. Shade Guide 000073)	1 Kit
000915	CC/Starter Kit	1 Kit
000920	CC/Introduction Kit	1 Kit
001590	CC/Gingiva Kit (inkl. Shade Guide 000048)	1 Kit
000250	CC/Bleach Shades Kit (inkl. Shade Guide 000049)	1 Kit
001503	CC/Shoulder Kit (inkl. Shade Guide 000047)	1 Kit
000801	CC/Kit 1 Transition Dentine (inkl. Shade Guide 000087)	1 Kit
620501	ZI-CT/Kit 1	1 Kit
620502	ZI-CT/Kit 2	1 Kit
620503	ZI-CT/Starter Kit	1 Kit
620506	ZI-CT/Introduction Kit	1 Kit
620504	ZI-CT/Gingiva Kit	1 Kit
620505	ZI-CT/Bleach Kit	1 Kit
620507	ZI-CT/Shoulder Kit	1 Kit
653010	ZI-F/Kit 1 (inkl. Shade Guide 000076)	1 Kit
653015	ZI-F/Kit 2 (inkl. Shade Guide 000077)	1 Kit
653025	ZI-F/Starter Kit	1 Kit
653030	ZI-F/Introduction Kit	1 Kit
603030	ZI-F/Gingiva Kit (inkl. Shade Guide 000048)	1 Kit
652050	ZI-F/Bleach Shades Kit (inkl. Shade Guide 000049)	1 Kit
830080	LS Press /Starter Kit	1 Kit
830081	LS Press /Introduction Kit A1, A2	1 Kit
830082	LS Press /Introduction Kit A3	1 Kit
830083	LS/Kit	1 Kit
001205	Magic Colour Professional Kit	1 Kit
001206	Magic Colour Starter Kit ZI-F	1 Kit
001207	Magic Colour Starter Kit LS	1 Kit
001200	CreaColor Kit	1 Kit
001100	Make up Neo Kit	1 Kit
001150	In Nova Neo Kit	1 Kit



Article no.	Article name	Content
	PASTES 3 g – 18 g	
000080	CC.LF/Crea Alloy Bond	2 x 4 g
000301	CC/Creapast A1	3 g
000302	CC/Creapast A2	3 g
000303	CC/Creapast A3	3 g
000304	CC/Creapast A3,5	3 g
000305	CC/Creapast A4	3 g
000306	CC/Creapast B1	3 g
000307	CC/Creapast B2	3 g
000308	CC/Creapast B3	3 g
000309	CC/Creapast B4	3 g
000310	CC/Creapast C1	3 g
000311	CC/Creapast C2	3 g
000312	CC/Creapast C3	3 g
000313	CC/Creapast C4	3 g
000314	CC/Creapast D2	3 g
000315	CC/Creapast D3	3 g
000316	CC/Creapast D4	3 g
000317	CC/Creapast IO white	3 g
000318	CC/Creapast IO violet	3 g
000319	CC/Creapast IO ochre	3 g
000320	CC/Creapast IO gold yellow	3 g
000322	CC/Creapast WOP sand	3 g
000323	CC/Creapast WOP sahara	3 g
000324	CC/Creapast WOP gold	3 g
000325	CC/Creapast G gingiva	3 g
000326	CC/Bleach Creapast O-AB	3 g
006301	CC/Creapast A1	18 g
006302	CC/Creapast A2	18 g
006303	CC/Creapast A3	18 g
006304	CC/Creapast A3,5	18 g
006305	CC/Creapast A4	18 g
006306	CC/Creapast B1	18 g
006307	CC/Creapast B2	18 g
006308	CC/Creapast B3	18 g
006309	CC/Creapast B4	18 g
006310	CC/Creapast C1	18 g
006311	CC/Creapast C2	18 g
006312	CC/Creapast C3	18 g
006313	CC/Creapast C4	18 g
006314	CC/Creapast D2	18 g
006315	CC/Creapast D3	18 g
006316	CC/Creapast D4	18 g
006317	CC/Creapast IO white	18 g
006318	CC/Creapast IO violet	18 g
006319	CC/Creapast IO ochre	18 g
006320	CC/Creapast IO gold yellow	18 g
006322	CC/Creapast WOP sand	18 g

Article no.	Article name	Content
006323	CC/Creapast WOP sahara	18 g
006324	CC/Creapast WOP gold	18 g
006325	CC/Creapast G gingiva	18 g
006326	CC/Bleach Creapast O-AB	18 g
	PORCELAIN 10 g – 15 g	
000174	CC/Glaze GL	10 g
001511	CC/Opaque Modifier OM-11 white	15 g
001512	CC/Opaque Modifier OM-12 pink	15 g
001513	CC/Opaque Modifier OM-13 lilac	15 g
001514	CC/Opaque Modifier OM-14 grey	15 g
001515	CC/Opaque Modifier OM-15 sahara	15 g
001516	CC/Opaque Modifier OM-16 ocker	15 g
001517	CC/Opaque Modifier OM-17 brown	15 g
001518	CC/Opaque Modifier OM-G gingiva	15 g
001521	CC/Shoulder SP-21 neutral	15 g
001522	CC/Shoulder SP-22 flamingo	15 g
001523	CC/Shoulder SP-23 sand	15 g
001524	CC/Shoulder SP-24 gold	15 g
001525	CC/Shoulder SP-25 red brown	15 g
001526	CC/Shoulder SP-26 yellow grey	15 g
001527	CC/Shoulder SP-27 ivory	15 g
001528	CC/Shoulder SP-28 opaque yellow	15 g
001529	CC/Shoulder SP-29 honey	15 g
001530	CC/Shoulder SP-G gingiva	15 g
000256	CC/Bleach Shoulder SP-AB	15 g
001532	CC/Opaque Dentine OD-32 (A2-A3)	15 g
001537	CC/Opaque Dentine OD-37 (B4)	15 g
001541	CC/Opaque Dentine OD-41 (D3)	15 g
001543	CC/Opaque Dentine OD-43 (ivory)	15 g
001544	CC/Opaque Dentine OD-44 (cuba)	15 g
001551	CC/Neck Transpa HT-51 iridescent	15 g
001552	CC/Neck Transpa HT-52 khaki	15 g
001553	CC/Neck Transpa HT-53 sand	15 g
001554	CC/Neck Transpa HT-54 honey	15 g
001555	CC/Neck Transpa HT-55 bordeaux	15 g
001556	CC/Neck Transpa HT-56 ochre	15 g
001561	CC/Make In MI-61 ivory	15 g
001562	CC/Make In MI-62 lemon	15 g
001563	CC/Make In MI-63 honey	15 g
001564	CC/Make In MI-64 flamingo	15 g
001565	CC/Make In MI-65 gold	15 g
001566	CC/Make In MI-66 olive	15 g
001568	CC/Approximal Dentine AD-1 light yellow	15 g
001569	CC/Approximal Dentine AD-2 orange	15 g
001570	CC/Opaque Dentine OD-A1	15 g
001571	CC/Opaque Dentine OD-A2	15 g
001572	CC/Opaque Dentine OD-A3	15 g

Article no.	Article name	Content
001573	CC/Opaque Dentine OD-A3,5	15 g
001574	CC/Opaque Dentine OD-A4	15 g
001575	CC/Opaque Dentine OD-B1	15 g
001576	CC/Opaque Dentine OD-B2	15 g
001577	CC/Opaque Dentine OD-B3	15 g
001578	CC/Opaque Dentine OD-B4	15 g
001579	CC/Opaque Dentine OD-C1	15 g
001580	CC/Opaque Dentine OD-C2	15 g
001581	CC/Opaque Dentine OD-C3	15 g
001582	CC/Opaque Dentine OD-C4	15 g
001583	CC/Opaque Dentine OD-D2	15 g
001584	CC/Opaque Dentine OD-D3	15 g
001585	CC/Opaque Dentine OD-D4	15 g
001591	CC/Gingiva G1 purple	15 g
001592	CC/Gingiva G2 dark pink	15 g
001593	CC/Gingiva G3 light pink	15 g
001594	CC/Gingiva G4 flamingo	15 g
001595	CC/Gingiva G5 rose	15 g
001596	CC/Gingiva G6 dark pink opaque	15 g
001597	CC/Gingiva GN neutral	15 g
001598	CC/Gingiva G7 violet	15 g
	PORCELAIN 20 g	
000201	CC/Opaque O-A1	20 g
000202	CC/Opaque O-A2	20 g
000203	CC/Opaque O-A3	20 g
000204	CC/Opaque O-A3,5	20 g
000205	CC/Opaque O-A4	20 g
000206	CC/Opaque O-B1	20 g
000207	CC/Opaque O-B2	20 g
000208	CC/Opaque O-B3	20 g
000209	CC/Opaque O-B4	20 g
000210	CC/Opaque O-C1	20 g
000211	CC/Opaque O-C2	20 g
000212	CC/Opaque O-C3	20 g
000213	CC/Opaque O-C4	20 g
000214	CC/Opaque O-D2	20 g
000215	CC/Opaque O-D3	20 g
000216	CC/Opaque O-D4	20 g
000217	CC/Opaque FOP	20 g
000231	CC/Dentine D-A1	20 g
000232	CC/Dentine D-A2	20 g
000233	CC/Dentine D-A3	20 g
000234	CC/Dentine D-A3,5	20 g
000235	CC/Dentine D-A4	20 g
000236	CC/Dentine D-B1	20 g
000237	CC/Dentine D-B2	20 g
000238	CC/Dentine D-B3	20 g

Article no.	Article name	Content
000239	CC/Dentine D-B4	20 g
000240	CC/Dentine D-C1	20 g
000241	CC/Dentine D-C2	20 g
000242	CC/Dentine D-C3	20 g
000243	CC/Dentine D-C4	20 g
000244	CC/Dentine D-D2	20 g
000245	CC/Dentine D-D3	20 g
000246	CC/Dentine D-D4	20 g
000247	CC/Dentine D-CW California white	20 g
000251	CC/Bleach Opaque O-AB	20 g
000252	CC/Bleach Dentine BD-A	20 g
000253	CC/Bleach Dentine BD-B	20 g
000254	CC/Bleach Dentine BD-B0	20 g
000255	CC/Bleach Enamel S-AB	20 g
000261	CC/Enamel S-057	20 g
000262	CC/Enamel S-058	20 g
000263	CC/Enamel S-059	20 g
000264	CC/Enamel S-060	20 g
000265	CC/Effect Enamel SI-01 light yellow	20 g
000266	CC/Effect Enamel SI-02 medium yellow	20 g
000267	CC/Effect Enamel SI-03 heavy yellow	20 g
000268	CC/Effect Enamel SI-04 light orange	20 g
000269	CC/Effect Enamel SI-05 medium orange	20 g
000270	CC/Effect Enamel SI-06 heavy orange	20 g
000271	CC/Window UC	20 g
000272	CC/Transpa NT neutral	20 g
000273	CC/Clear CL-O	20 g
000274	CC/Glaze GL	20 g
000287	CC/Transpa OT+	20 g
000288	CC/Transpa Geller-Mix	20 g
000275	CC/Transpa OT opal	20 g
000276	CC/Transpa TI-1 blue	20 g
000277	CC/Transpa TI-2 white	20 g
000278	CC/Transpa TI-3 pink	20 g
000279	CC/Transpa TI-4 yellow	20 g
000280	CC/Transpa TI-5 grey	20 g
000281	CC/Pearl Enamel PS-0 white	20 g
000282	CC/Pearl Enamel PS-1 blue/white	20 g
000283	CC/Pearl Enamel PS-2 flamingo	20 g
000284	CC/Pearl Enamel PS-3 orange	20 g
000285	CC/Effect Enamel SO-10 blue	20 g
000286	CC/Effect Enamel SO-11 orange	20 g
000290	CC/Correction Powder KM	20 g
003201	CC/Transition Dentine TD-A1	20 g
003202	CC/Transition Dentine TD-A2	20 g
003203	CC/Transition Dentine TD-A3	20 g
003204	CC/Transition Dentine TD-A3,5	20 g
003205	CC/Transition Dentine TD-A4	20 g

new

Article no.	Article name	Content
003206	CC/Transition Dentine TD-B1	20 g
003207	CC/Transition Dentine TD-B2	20 g
003208	CC/Transition Dentine TD-B3	20 g
003209	CC/Transition Dentine TD-B4	20 g
003210	CC/Transition Dentine TD-C1	20 g
003211	CC/Transition Dentine TD-C2	20 g
003212	CC/Transition Dentine TD-C3	20 g
003213	CC/Transition Dentine TD-C4	20 g
003214	CC/Transition Dentine TD-D2	20 g
003215	CC/Transition Dentine TD-D3	20 g
003216	CC/Transition Dentine TD-D4	20 g
003217	CC/Transition Dentine TD-BA	20 g
	PORCELAIN 50 g	
000501	CC/Opaque O-A1	50 g
000502	CC/Opaque O-A2	50 g
000503	CC/Opaque O-A3	50 g
000504	CC/Opaque O-A3,5	50 g
000505	CC/Opaque O-A4	50 g
000506	CC/Opaque O-B1	50 g
000507	CC/Opaque O-B2	50 g
000508	CC/Opaque O-B3	50 g
000509	CC/Opaque O-B4	50 g
000510	CC/Opaque O-C1	50 g
000511	CC/Opaque O-C2	50 g
000512	CC/Opaque O-C3	50 g
000513	CC/Opaque O-C4	50 g
000514	CC/Opaque O-D2	50 g
000515	CC/Opaque O-D3	50 g
000516	CC/Opaque O-D4	50 g
000517	CC/Opaque O-FOP	50 g
000531	CC/Dentine D-A1	50 g
000532	CC/Dentine D-A2	50 g
000533	CC/Dentine D-A3	50 g
000534	CC/Dentine D-A3,5	50 g
000535	CC/Dentine D-A4	50 g
000536	CC/Dentine D-B1	50 g
000537	CC/Dentine D-B2	50 g
000538	CC/Dentine D-B3	50 g
000539	CC/Dentine D-B4	50 g
000540	CC/Dentine D-C1	50 g
000541	CC/Dentine D-C2	50 g
000542	CC/Dentine D-C3	50 g
000543	CC/Dentine D-C4	50 g
000544	CC/Dentine D-D2	50 g
000545	CC/Dentine D-D3	50 g
000546	CC/Dentine D-D4	50 g
000547	CC/Dentine D-CW California white	50 g

Article no.	Article name	Content
000551	CC/Bleach Opaque O-AB	50 g
000552	CC/Bleach Dentine BD-A	50 g
000553	CC/Bleach Dentine BD-B	50 g
000554	CC/Bleach Dentine BD-B0	50 g
000555	CC/Bleach Enamel S-AB	50 g
000556	CC/Bleach Shoulder SP-AB	50 g
000561	CC/Enamel S-057	50 g
000562	CC/Enamel S-058	50 g
000563	CC/Enamel S-059	50 g
000564	CC/Enamel S-060	50 g
000565	CC/Effect Enamel SI-01 light yellow	50 g
000566	CC/Effect Enamel SI-02 medium yellow	50 g
000567	CC/Effect Enamel SI-03 heavy yellow	50 g
000568	CC/Effect Enamel SI-04 light orange	50 g
000569	CC/Effect Enamel SI-05 medium orange	50 g
000570	CC/Effect Enamel SI-06 heavy orange	50 g
000571	CC/Window UC	50 g
000572	CC/Transpa NT neutral	50 g
000573	CC/Clear CL-O	50 g
000587	CC/Transpa OT+	50 g
000588	CC/Transpa Geller-Mix	50 g
000575	CC/Transpa OT opal	50 g
000576	CC/Transpa TI-1 blue	50 g
000577	CC/Transpa TI-2 white	50 g
000578	CC/Transpa TI-3 pink	50 g
000579	CC/Transpa TI-4 yellow	50 g
000580	CC/Transpa TI-5 grey	50 g
000581	CC/Pearl Enamel PS-0 white	50 g
000582	CC/Pearl Enamel PS-1 blue/white	50 g
000583	CC/Pearl Enamel PS-2 flamingo	50 g
000584	CC/Pearl Enamel PS-3 orange	50 g
000585	CC/Effect Enamel SO-10 blue	50 g
000586	CC/Effect Enamel SO-11 orange	50 g
003501	CC/Transition Dentine TD-A1	50 g
003502	CC/Transition Dentine TD-A2	50 g
003503	CC/Transition Dentine TD-A3	50 g
003504	CC/Transition Dentine TD-A3,5	50 g
003505	CC/Transition Dentine TD-A4	50 g
003506	CC/Transition Dentine TD-B1	50 g
003507	CC/Transition Dentine TD-B2	50 g
003508	CC/Transition Dentine TD-B3	50 g
003509	CC/Transition Dentine TD-B4	50 g
003510	CC/Transition Dentine TD-C1	50 g
003511	CC/Transition Dentine TD-C2	50 g
003512	CC/Transition Dentine TD-C3	50 g
003513	CC/Transition Dentine TD-C4	50 g
003514	CC/Transition Dentine TD-D2	50 g
003515	CC/Transition Dentine TD-D3	50 g

new

Article no.	Article name	Content
003516	CC/Transition Dentine TD-D4	50 g
003517	CC/Transition Dentine TD-BA	50 g
005011	CC/Opaque Modifier OM-11 white	50 g
005012	CC/Opaque Modifier OM-12 pink	50 g
005013	CC/Opaque Modifier OM-13 lilac	50 g
005014	CC/Opaque Modifier OM-14 grey	50 g
005015	CC/Opaque Modifier OM-15 sahara	50 g
005016	CC/Opaque Modifier OM-16 ocker	50 g
005017	CC/Opaque Modifier OM-17 brown	50 g
005018	CC/Opaque Modifier OM-G gingiva	50 g
005021	CC/Shoulder SP-21 neutral	50 g
005022	CC/Shoulder SP-22 flamingo	50 g
005023	CC/Shoulder SP-23 sand	50 g
005024	CC/Shoulder SP-24 gold	50 g
005025	CC/Shoulder SP-25 red brown	50 g
005026	CC/Shoulder SP-26 yellow grey	50 g
005027	CC/Shoulder SP-27 ivory	50 g
005028	CC/Shoulder SP-28 opaque yellow	50 g
005029	CC/Shoulder SP-29 honey	50 g
005030	CC/Shoulder SP-G gingiva	50 g
005032	CC /Opaque Dentine OD-32 (A2-A3)	50 g
005037	CC /Opaque Dentine OD-37 (B4)	50 g
005041	CC /Opaque Dentine OD-41 (D3)	50 g
005043	CC /Opaque Dentine OD-43 (ivory)	50 g
005044	CC /Opaque Dentine OD-44 (cuba)	50 g
005051	CC /Neck Transpa HT-51 iridescent	50 g
005052	CC/Neck Transpa HT-52 khaki	50 g
005053	CC/Neck Transpa HT-53 sand	50 g
005054	CC/Neck Transpa HT-54 honey	50 g
005055	CC/Neck Transpa HT-55 bordeaux	50 g
005056	CC/Neck Transpa HT-56 ochre	50 g
005061	CC/Make In MI-61 ivory	50 g
005062	CC/Make In MI-62 lemon	50 g
005063	CC/Make In MI-63 honey	50 g
005064	CC/Make In MI-64 flamingo	50 g
005065	CC/Make In MI-65 gold	50 g
005066	CC/Make In MI-66 olive	50 g
005068	CC/Approximal Dentine AD-1 light yellow	50 g
005069	CC/Approximal Dentine AD-2 orange	50 g
005070	CC/Opaque Dentine OD-A1	50 g
005071	CC/Opaque Dentine OD-A2	50 g
005072	CC/Opaque Dentine OD-A3	50 g
005073	CC/Opaque Dentine OD-A3,5	50 g
005074	CC/Opaque Dentine OD-A4	50 g
005075	CC/Opaque Dentine OD-B1	50 g
005076	CC/Opaque Dentine OD-B2	50 g
005077	CC/Opaque Dentine OD-B3	50 g
005078	CC/Opaque Dentine OD-B4	50 g

Article no.	Article name	Content
005079	CC/Opaque Dentine OD-C1	50 g
005080	CC/Opaque Dentine OD-C2	50 g
005081	CC/Opaque Dentine OD-C3	50 g
005082	CC/Opaque Dentine OD-C4	50 g
005083	CC/Opaque Dentine OD-D2	50 g
005084	CC/Opaque Dentine OD-D3	50 g
005085	CC/Opaque Dentine OD-D4	50 g
005591	CC/Gingiva G1 purple	50 g
005592	CC/Gingiva G2 dark pink	50 g
005593	CC/Gingiva G3 light pink	50 g
005594	CC/Gingiva G4 flamingo	50 g
005595	CC/Gingiva G5 rose	50 g
005596	CC/Gingiva G6 dark pink opaque	50 g
005597	CC/Gingiva GN neutral	50 g
005598	CC/Gingiva G7 violet	50 g
PORCELAIN 250 g		
002501	CC/Opaque O-A1	250 g
002502	CC/Opaque O-A2	250 g
002503	CC/Opaque O-A3	250 g
002504	CC/Opaque O-A3,5	250 g
002505	CC/Opaque O-A4	250 g
002506	CC/Opaque O-B1	250 g
002507	CC/Opaque O-B2	250 g
002508	CC/Opaque O-B3	250 g
002509	CC/Opaque O-B4	250 g
002510	CC/Opaque O-C1	250 g
002511	CC/Opaque O-C2	250 g
002512	CC/Opaque O-C3	250 g
002513	CC/Opaque O-C4	250 g
002514	CC/Opaque O-D2	250 g
002515	CC/Opaque O-D3	250 g
002516	CC/Opaque O-D4	250 g
002531	CC/Dentine D-A1	250 g
002532	CC/Dentine D-A2	250 g
002533	CC/Dentine D-A3	250 g
002534	CC/Dentine D-A3,5	250 g
002535	CC/Dentine D-A4	250 g
002536	CC/Dentine D-B1	250 g
002537	CC/Dentine D-B2	250 g
002538	CC/Dentine D-B3	250 g
002539	CC/Dentine D-B4	250 g
002540	CC/Dentine D-C1	250 g
002541	CC/Dentine D-C2	250 g
002542	CC/Dentine D-C3	250 g
002543	CC/Dentine D-C4	250 g
002544	CC/Dentine D-D2	250 g
002545	CC/Dentine D-D3	250 g

Article no.	Article name	Content
002546	CC/Dentine D-D4	250 g
002547	CC/California white	250 g
002561	CC/Enamel S-057	250 g
002562	CC/Enamel S-058	250 g
002563	CC/Enamel S-059	250 g
002564	CC/Enamel S-060	250 g
002565	CC/Effect Enamel SI-01 light yellow	250 g
002566	CC/Effect Enamel SI-02 medium yellow	250 g
002567	CC/Effect Enamel SI-03 heavy yellow	250 g
002568	CC/Effect Enamel SI-04 light orange	250 g
002569	CC/Effect Enamel SI-05 medium orange	250 g
002570	CC/Effect Enamel SI-06 heavy orange	250 g
002572	CC/Transpa NT neutral	250 g
002573	CC/Clear CL-O	250 g
002575	CC/Transpa OT opal	250 g
002576	CC/Transpa TI-1 blue	250 g
002577	CC/Transpa TI-2 white	250 g
002578	CC/Transpa TI-3 pink	250 g
002579	CC/Transpa TI-4 yellow	250 g
002580	CC/Transpa TI-5 grey	250 g
002581	CC/Pearl Enamel PS-0 white	250 g
002582	CC/Pearl Enamel PS-1 blue/white	250 g
002583	CC/Pearl Enamel PS-2 flamingo	250 g
002584	CC/Pearl Enamel PS-3 orange	250 g
025071	CC/Opaque Dentine OD-A2	250 g
025072	CC/Opaque Dentine OD-A3	250 g
025073	CC/Opaque Dentine OD-A3,5	250 g
025076	CC/Opaque Dentine OD-B2	250 g
025077	CC/Opaque Dentine OD-B3	250 g
LIQUIDS		
000012	Modelling Liquid high-fusing ML-H	25 ml
000013	Modelling Liquid high-fusing ML-H	50 ml
000014	Modelling Liquid high-fusing ML-H	250 ml
000015	Modelling Liquid high-fusing ML-H	1000 ml
000019	Opaque Liquid high-fusing OL-H	25 ml
000020	Opaque Liquid high-fusing OL-H	50 ml
000021	Opaque Liquid high-fusing OL-H	250 ml
000010	Universal & Glaze Liquid UF/GL	25 ml
000011	Universal & Glaze Liquid UF/GL	50 ml
000025	Shoulder Liquid SPL	25 ml
000026	Shoulder Liquid SPL	50 ml
000090	Moisturing Fluid	25 ml

CREATION ZI-CT | ZIRCONIUM DIOXIDE CERAMIC

Article no.	Article name	Content
	PORCELAIN 4 g – 15 g	
620061	ZI-CT/Frame Shade NT	4 g
620074	ZI-CT/Glaze GL	10 g
620101	ZI-CT/Gingiva G1 purple	15 g
620102	ZI-CT/Gingiva G2 dark pink	15 g
620103	ZI-CT/Gingiva G3 light pink	15 g
620104	ZI-CT/Gingiva G4 flamingo	15 g
620105	ZI-CT/Gingiva G5 rose	15 g
620106	ZI-CT/Gingiva G6 dark pink opaque	15 g
620107	ZI-CT/Gingiva G7 violet	15 g
620108	ZI-CT/Gingiva GN neutral	15 g
620120	ZI-CT/Shoulder SP-AB	15 g
620121	ZI-CT/Shoulder SP-21 neutral	15 g
620122	ZI-CT/Shoulder SP-22 flamingo	15 g
620123	ZI-CT/Shoulder SP-23 sand	15 g
620124	ZI-CT/Shoulder SP-24 gold	15 g
620125	ZI-CT/Shoulder SP-25 red brown	15 g
620126	ZI-CT/Shoulder SP-26 yellow grey	15 g
620127	ZI-CT/Shoulder SP-27 ivory	15 g
620128	ZI-CT/Shoulder SP-28 opaque yellow	15 g
620129	ZI-CT/Shoulder SP-29 honey	15 g
620131	ZI-CT/Opaque Dentine OD-32 havanna	15 g
620132	ZI-CT/Opaque Dentine OD-37 curry	15 g
620133	ZI-CT/Opaque Dentine OD-41 orange	15 g
620134	ZI-CT/Opaque Dentine OD-43 ivory	15 g
620135	ZI-CT/Opaque Dentine OD-44 cuba	15 g
620146	ZI-CT/Approximal Dentine AD-1 light yellow	15 g
620147	ZI-CT/Approximal Dentine AD-2 orange	15 g
620151	ZI-CT/Neck Transpa HT-51 iridescent	15 g
620152	ZI-CT/Neck Transpa HT-52 kaki	15 g
620153	ZI-CT/Neck Transpa HT-53 sand	15 g
620154	ZI-CT/Neck Transpa HT-54 honey	15 g
620155	ZI-CT/Neck Transpa HT-55 bordeaux	15 g
620156	ZI-CT/Neck Transpa HT-56 ochre	15 g
620161	ZI-CT/Make In MI-61 ivory	15 g
620162	ZI-CT/Make In MI-62 lemon	15 g
620163	ZI-CT/Make In MI-63 honey	15 g
620164	ZI-CT/Make In MI-64 flamingo	15 g
620165	ZI-CT/Make In MI-65 gold	15 g
620166	ZI-CT/Make In MI-66 olive	15 g
	PORCELAIN 20 g	
620001	ZI-CT/Dentine D-A1	20 g
620002	ZI-CT/Dentine D-A2	20 g
620003	ZI-CT/Dentine D-A3	20 g
620004	ZI-CT/Dentine D-A3,5	20 g
620005	ZI-CT/Dentine D-A4	20 g
620006	ZI-CT/Dentine D-B1	20 g

CREATION ZI-CT | ZIRCONIUM DIOXIDE CERAMIC

Article no.	Article name	Content
620007	ZI-CT/Dentine D-B2	20 g
620008	ZI-CT/Dentine D-B3	20 g
620009	ZI-CT/Dentine D-B4	20 g
620010	ZI-CT/Dentine D-C1	20 g
620011	ZI-CT/Dentine D-C2	20 g
620012	ZI-CT/Dentine D-C3	20 g
620013	ZI-CT/Dentine D-C4	20 g
620014	ZI-CT/Dentine D-D2	20 g
620015	ZI-CT/Dentine D-D3	20 g
620016	ZI-CT/Dentine D-D4	20 g
620017	ZI-CT/Enamel E-57	20 g
620018	ZI-CT/Enamel E-58	20 g
620019	ZI-CT/Enamel E-59	20 g
620020	ZI-CT/Enamel E-60	20 g
620021	ZI-CT/Clear CL-O	20 g
620022	ZI-CT/Window UC	20 g
620023	ZI-CT/Bleach Dentine BD-A	20 g
620024	ZI-CT/Bleach Dentine BD-B	20 g
620025	ZI-CT/Bleach Dentine BD-B0	20 g
620026	ZI-CT/Bleach Enamel S-AB	20 g
620031	ZI-CT/Effect Enamel SI-01 light yellow	20 g
620032	ZI-CT/Effect Enamel SI-02 medium yellow	20 g
620033	ZI-CT/Effect Enamel SI-03 heavy yellow	20 g
620034	ZI-CT/Effect Enamel SI-04 light orange	20 g
620035	ZI-CT/Effect Enamel SI-05 medium orange	20 g
620036	ZI-CT/Effect Enamel SI-06 heavy orange	20 g
620037	ZI-CT/Effect Enamel SO-10 blue	20 g
620038	ZI-CT/Effect Enamel SO-11 orange	20 g
620039	ZI-CT/Transpa NT neutral	20 g
620046	ZI-CT/Transpa OT+	20 g
620047	ZI-CT/Transpa Geller-Mix	20 g
620040	ZI-CT/Transpa OT opal	20 g
620041	ZI-CT/Transpa TI-1 blue	20 g
620042	ZI-CT/Transpa TI-2 white	20 g
620043	ZI-CT/Transpa TI-3 pink	20 g
620044	ZI-CT/Transpa TI-4 yellow	20 g
620045	ZI-CT/Transpa TI-5 grey	20 g
620050	ZI-CT/Pearl Enamel PS-0 white	20 g
620051	ZI-CT/Pearl Enamel PS-1 blue/white	20 g
620052	ZI-CT/Pearl Enamel PS-2 flamingo	20 g
620053	ZI-CT/Pearl Enamel PS-3 orange	20 g
620060	ZI-CT/Correction Powder KM	20 g
	PORCELAIN 50 g	
620201	ZI-CT/Dentine D-A1	50 g
620202	ZI-CT/Dentine D-A2	50 g
620203	ZI-CT/Dentine D-A3	50 g
620204	ZI-CT/Dentine D-A3,5	50 g



CREATION ZI-CT | ZIRCONIUM DIOXIDE CERAMIC

Article no.	Article name	Content
620205	ZI-CT/Dentine D-A4	50 g
620206	ZI-CT/Dentine D-B1	50 g
620207	ZI-CT/Dentine D-B2	50 g
620208	ZI-CT/Dentine D-B3	50 g
620209	ZI-CT/Dentine D-B4	50 g
620210	ZI-CT/Dentine D-C1	50 g
620211	ZI-CT/Dentine D-C2	50 g
620212	ZI-CT/Dentine D-C3	50 g
620213	ZI-CT/Dentine D-C4	50 g
620214	ZI-CT/Dentine D-D2	50 g
620215	ZI-CT/Dentine D-D3	50 g
620216	ZI-CT/Dentine D-D4	50 g
620217	ZI-CT/Enamel E-57	50 g
620218	ZI-CT/Enamel E-58	50 g
620219	ZI-CT/Enamel E-59	50 g
620220	ZI-CT/Enamel E-60	50 g
620221	ZI-CT/Clear CL-O	50 g
620222	ZI-CT/Window UC	50 g
620223	ZI-CT/Bleach Dentine BD-A	50 g
620224	ZI-CT/Bleach Dentine BD-B	50 g
620225	ZI-CT/Bleach Dentine BD-B0	50 g
620226	ZI-CT/Bleach Enamel S-AB	50 g
620231	ZI-CT/Effect Enamel SI-01 light yellow	50 g
620232	ZI-CT/Effect Enamel SI-02 medium yellow	50 g
620233	ZI-CT/Effect Enamel SI-03 heavy yellow	50 g
620234	ZI-CT/Effect Enamel SI-04 light orange	50 g
620235	ZI-CT/Effect Enamel SI-05 medium orange	50 g
620236	ZI-CT/Effect Enamel SI-06 heavy orange	50 g
620237	ZI-CT/Effect Enamel SO-10 blue	50 g
620238	ZI-CT/Effect Enamel SO-11 orange	50 g
620239	ZI-CT/Transpa NT neutral	50 g
620246	ZI-CT/Transpa OT+	50 g
620247	ZI-CT/Transpa Geller-Mix	50 g
620240	ZI-CT/Transpa OT opal	50 g
620241	ZI-CT/Transpa TI-1 blue	50 g
620242	ZI-CT/Transpa TI-2 white	50 g
620243	ZI-CT/Transpa TI-3 pink	50 g
620244	ZI-CT/Transpa TI-4 yellow	50 g
620245	ZI-CT/Transpa TI-5 grey	50 g
620250	ZI-CT/Pearl Enamel PS-0 white	50 g
620251	ZI-CT/Pearl Enamel PS-1 blue/white	50 g
620252	ZI-CT/Pearl Enamel PS-2 flamingo	50 g
620253	ZI-CT/Pearl Enamel PS-3 orange	50 g
620261	ZI-CT/Opaque Dentine OD-32 havanna	50 g
620262	ZI-CT/Opaque Dentine OD-37 curry	50 g
620263	ZI-CT/Opaque Dentine OD-41 orange	50 g
620264	ZI-CT/Opaque Dentine OD-43 ivory	50 g
620265	ZI-CT/Opaque Dentine OD-44 cuba	50 g



CREATION ZI-CT | ZIRCONIUM DIOXIDE CERAMIC

Article no.	Article name	Content
620266	ZI-CT/Approximal Dentine AD-1 light yellow	50 g
620267	ZI-CT/Approximal Dentine AD-2 orange	50 g
620271	ZI-CT/Neck Transpa HT-51 iridescent	50 g
620272	ZI-CT/Neck Transpa HT-52 kaki	50 g
620273	ZI-CT/Neck Transpa HT-53 sand	50 g
620274	ZI-CT/Neck Transpa HT-54 honey	50 g
620275	ZI-CT/Neck Transpa HT-55 bordeaux	50 g
620276	ZI-CT/Neck Transpa HT-56 ochre	50 g
	PORCELAIN 250 g	
620302	ZI-CT/Dentine D-A2	250 g
620303	ZI-CT/Dentine D-A3	250 g
620318	ZI-CT/Enamel E-58	250 g
620319	ZI-CT/Enamel E-59	250 g
620321	ZI-CT/Clear CL-O	250 g
	LIQUIDS	
000012	Modelling Liquid high-fusing ML-H	25 ml
000013	Modelling Liquid high-fusing ML-H	50 ml
000014	Modelling Liquid high-fusing ML-H	250 ml
000015	Modelling Liquid high-fusing ML-H	1000 ml
000010	Universal & Glaze Liquid UF/GL	25 ml
000011	Universal & Glaze Liquid UF/GL	50 ml
000025	Shoulder Liquid SPL	25 ml
000026	Shoulder Liquid SPL	50 ml
000090	Moisturing Fluid	25 ml

CREATION ZI-F | LOW-FUSING ZIRCONIUM DIOXIDE CERAMIC

Article no.	Article name	Content
	PORCELAIN 4g – 15 g	
601087	ZI-F/Frame Shade NT	4 g
601074	ZI-F/Glaze GL	10 g
601081	ZI-F/Frame Shade FS-1	10 g
601082	ZI-F/Frame Shade FS-2	10 g
601083	ZI-F/Frame Shade FS-3	10 g
601084	ZI-F/Frame Shade FS-4	10 g
601085	ZI-F/Frame Shade FS-IL (Illusion)	10 g
601086	ZI-F/Frame Shade FS-G (Gingiva)	10 g
602051	ZI-F/Bleach Frame Shade FS-AB	10 g
601091	ZI-F/Neck Dentine HD-01	15 g
601092	ZI-F/Neck Dentine HD-02	15 g
601093	ZI-F/Neck Dentine HD-03	15 g
601094	ZI-F/Neck Dentine HD-04	15 g
602056	ZI-F/Bleach Shoulder SP-AB	15 g
602121	ZI-F/Shoulder SP-21 neutral	15 g
602122	ZI-F/Shoulder SP-22 flamingo	15 g
602123	ZI-F/Shoulder SP-23 sand	15 g
602124	ZI-F/Shoulder SP-24 gold	15 g
602125	ZI-F/Shoulder SP-25 red brown	15 g
602128	ZI-F/Shoulder SP-28 opaque yellow	15 g
602145	ZI-F/Gingiva G1 purple	15 g
602146	ZI-F/Gingiva G2 dark pink	15 g
602147	ZI-F/Gingiva G3 light pink	15 g
602148	ZI-F/Gingiva G4 flamingo	15 g
602149	ZI-F/Gingiva G5 rose	15 g
602150	ZI-F/Gingiva G6 dark pink opaque	15 g
602151	ZI-F/Gingiva GN neutral	15 g
602157	ZI-F/Gingiva G7 violet	15 g
602161	ZI-F/Make In MI-61 ivory	15 g
602163	ZI-F/Make In MI-63 honey	15 g
602165	ZI-F/Make In MI-65 gold	15 g
652132	ZI-F/Opaque Dentine OD-32 havanna	15 g
652137	ZI-F/Opaque Dentine OD-37 curry	15 g
652141	ZI-F/Opaque Dentine OD-41 orange	15 g
652143	ZI-F/Opaque Dentine OD-43 ivory	15 g
652144	ZI-F/Opaque Dentine OD-44 cuba	15 g
652151	ZI-F/Neck Transpa HT-51 iridescent	15 g
652152	ZI-F/Neck Transpa HT-52 kaki	15 g
652153	ZI-F/Neck Transpa HT-53 sand	15 g
652154	ZI-F/Neck Transpa HT-54 honey	15 g
652156	ZI-F/Neck Transpa HT-56 ochre	15 g
652168	ZI-F/Approximal Dentine AD-1 light yellow	15 g
652169	ZI-F/Approximal Dentine AD-2 orange	15 g
	PORCELAIN 20 g	
652001	ZI-F/Dentine D-A1	20 g
652002	ZI-F/Dentine D-A2	20 g

CREATION ZI-F | LOW-FUSING ZIRCONIUM DIOXIDE CERAMIC

Article no.	Article name	Content
652003	ZI-F/Dentine D-A3	20 g
652004	ZI-F/Dentine D-A3,5	20 g
652005	ZI-F/Dentine D-A4	20 g
652006	ZI-F/Dentine D-B1	20 g
652007	ZI-F/Dentine D-B2	20 g
652008	ZI-F/Dentine D-B3	20 g
652009	ZI-F/Dentine D-B4	20 g
652010	ZI-F/Dentine D-C1	20 g
652011	ZI-F/Dentine D-C2	20 g
652012	ZI-F/Dentine D-C3	20 g
652013	ZI-F/Dentine D-C4	20 g
652014	ZI-F/Dentine D-D2	20 g
652015	ZI-F/Dentine D-D3	20 g
652016	ZI-F/Dentine D-D4	20 g
652052	ZI-F/Bleach Dentine BD-A	20 g
652053	ZI-F/Bleach Dentine BD-B	20 g
652054	ZI-F/Bleach Dentine BD-B0	20 g
652055	ZI-F/Bleach Enamel S-AB	20 g
652057	ZI-F/Enamel E-57	20 g
652058	ZI-F/Enamel E-58	20 g
652059	ZI-F/Enamel E-59	20 g
652060	ZI-F/Enamel E-60	20 g
652066	ZI-F/Effect Enamel SI-02 medium yellow	20 g
652068	ZI-F/Effect Enamel SI-04 light orange	20 g
652070	ZI-F/Effect Enamel SI-06 heavy orange	20 g
602072	ZI-F/Window UC	20 g
652073	ZI-F/Clear CL-O	20 g
652074	ZI-F/Transpa NT neutral	20 g
652087	ZI-F/Transpa OT+	20 g
652088	ZI-F/Transpa Geller-Mix	20 g
652075	ZI-F/Transpa OT opal	20 g
652076	ZI-F/Transpa TI-1 blue	20 g
652077	ZI-F/Transpa TI-2 white	20 g
652078	ZI-F/Transpa TI-3 pink	20 g
652079	ZI-F/Transpa TI-4 yellow	20 g
652080	ZI-F/Transpa TI-5 grey	20 g
652081	ZI-F/Effect Enamel PS-0 white	20 g
652083	ZI-F/Effect Enamel PS-3 orange	20 g
652085	ZI-F/Effect Enamel SO-10 blue	20 g
652086	ZI-F/Effect Enamel SO-11 orange	20 g
652090	ZI-F/Correction Powder KM	20 g
602174	ZI-F/Glaze GL	20 g
	PORCELAIN 50 g	
605056	ZI-F/Bleach Shoulder SP-AB	50 g
605121	ZI-F/Shoulder SP-21 neutral	50 g
605122	ZI-F/Shoulder SP-22 flamingo	50 g
605123	ZI-F/Shoulder SP-23 sand	50 g



CREATION ZI-F | LOW-FUSING ZIRCONIUM DIOXIDE CERAMIC

Article no.	Article name	Content
605124	ZI-F/Shoulder SP-24 gold	50 g
605125	ZI-F/Shoulder SP-25 red brown	50 g
605128	ZI-F/Shoulder SP-28 opaque yellow	50 g
605145	ZI-F/Gingiva G1 purple	50 g
605146	ZI-F/Gingiva G2 dark pink	50 g
605147	ZI-F/Gingiva G3 light pink	50 g
605148	ZI-F/Gingiva G4 flamingo	50 g
605149	ZI-F/Gingiva G5 rose	50 g
605150	ZI-F/Gingiva G6 dark pink opaque	50 g
605151	ZI-F/Gingiva GN neutral	50 g
605157	ZI-F/Gingiva G7 violet	50 g
605161	ZI-F/Make In MI-61 ivory	50 g
605163	ZI-F/Make In MI-63 honey	50 g
605165	ZI-F/Make In MI-65 gold	50 g
655001	ZI-F/Dentine D-A1	50 g
655002	ZI-F/Dentine D-A2	50 g
655003	ZI-F/Dentine D-A3	50 g
655004	ZI-F/Dentine D-A3,5	50 g
655005	ZI-F/Dentine D-A4	50 g
655006	ZI-F/Dentine D-B1	50 g
655007	ZI-F/Dentine D-B2	50 g
655008	ZI-F/Dentine D-B3	50 g
655009	ZI-F/Dentine D-B4	50 g
655010	ZI-F/Dentine D-C1	50 g
655011	ZI-F/Dentine D-C2	50 g
655012	ZI-F/Dentine D-C3	50 g
655013	ZI-F/Dentine D-C4	50 g
655014	ZI-F/Dentine D-D2	50 g
655015	ZI-F/Dentine D-D3	50 g
655016	ZI-F/Dentine D-D4	50 g
655052	ZI-F/Bleach Dentine BD-A	50 g
655053	ZI-F/Bleach Dentine BD-B	50 g
655054	ZI-F/Bleach Dentine BD-B0	50 g
655055	ZI-F/Bleach Enamel S-AB	50 g
655057	ZI-F/Enamel E-57	50 g
655058	ZI-F/Enamel E-58	50 g
655059	ZI-F/Enamel E-59	50 g
655060	ZI-F/Enamel E-60	50 g
655066	ZI-F/Effect Enamel SI-02 medium yellow	50 g
655068	ZI-F/Effect Enamel SI-04 light orange	50 g
655070	ZI-F/Effect Enamel SI-06 heavy orange	50 g
605072	ZI-F/Window UC	50 g
655073	ZI-F/Clear CL-O	50 g
655074	ZI-F/Transpa NT neutral	50 g
655087	ZI-F/Transpa OT+	50 g
655088	ZI-F/Transpa Geller-Mix	50 g
655075	ZI-F/Transpa OT opal	50 g
655076	ZI-F/Transpa TI-1 blue	50 g



CREATION ZI-F | LOW-FUSING ZIRCONIUM DIOXIDE CERAMIC

Article no.	Article name	Content
655077	ZI-F/Transpa TI-2 white	50 g
655078	ZI-F/Transpa TI-3 pink	50 g
655079	ZI-F/Transpa TI-4 yellow	50 g
655080	ZI-F/Transpa TI-5 grey	50 g
655081	ZI-F/Effect Enamel PS-0 white	50 g
655083	ZI-F/Effect Enamel PS-3 orange	50 g
655085	ZI-F/Effect Enamel SO-10 blue	50 g
655086	ZI-F/Effect Enamel SO-11 orange	50 g
655132	ZI-F/Opaque Dentine OD-32 havanna	50 g
655137	ZI-F/Opaque Dentine OD-37 curry	50 g
655141	ZI-F/Opaque Dentine OD-41 orange	50 g
655143	ZI-F/Opaque Dentine OD-43 ivory	50 g
655144	ZI-F/Opaque Dentine OD-44 cuba	50 g
655151	ZI-F/Neck Transpa HT-51 iridescent	50 g
655152	ZI-F/Neck Transpa HT-52 kaki	50 g
655153	ZI-F/Neck Transpa HT-53 sand	50 g
655154	ZI-F/Neck Transpa HT-54 honey	50 g
655156	ZI-F/Neck Transpa HT-56 ochre	50 g
655168	ZI-F/Approximal Dentine AD-1 light yellow	50 g
655169	ZI-F/Approximal Dentine AD-2 orange	50 g
LIQUIDS		
000016	Modelling Liquid low-fusing ML-L	25 ml
000017	Modelling Liquid low-fusing ML-L	50 ml
000018	Modelling Liquid low-fusing ML-L	250 ml
000010	Universal & Glaze Liquid UF/GL	25 ml
000011	Universal & Glaze Liquid UF/GL	50 ml
000025	Shoulder Liquid SPL	25 ml
000026	Shoulder Liquid SPL	50 ml
000090	Moisturing Fluid	25 ml



Article no.	Article name	Content
	INGOTS 3 g	
830101	LS/LS Press, MO-0	5 x 3 g
830102	LS/LS Press, MO-1	5 x 3 g
830103	LS/LS Press, MO-2	5 x 3 g
830129	LS/LS Press LT-B00	5 x 3 g
830128	LS/LS Press LT-B0	5 x 3 g
830124	LS/LS Press LT-A1	5 x 3 g
830125	LS/LS Press LT-A2	5 x 3 g
830126	LS/LS Press LT-A3	5 x 3 g
830127	LS/LS Press LT-B1	5 x 3 g
830104	LS/LS Press LT-A	5 x 3 g
830105	LS/LS Press LT-B	5 x 3 g
830106	LS/LS Press LT-C	5 x 3 g
830107	LS/LS Press LT-D	5 x 3 g
830108	LS/LS Press MT-B00	5 x 3 g
830109	LS/LS Press MT-B0	5 x 3 g
830110	LS/LS Press MT-A1	5 x 3 g
830111	LS/LS Press MT-A2	5 x 3 g
830112	LS/LS Press MT-A3	5 x 3 g
830113	LS/LS Press MT-B1	5 x 3 g
830114	LS/LS Press MT-B2	5 x 3 g
830117	LS/LS Press MT-D2	5 x 3 g
830122	LS/LS Press HT-B00	5 x 3 g
830123	LS/LS Press HT-B0	5 x 3 g
830118	LS/LS Press HT-1	5 x 3 g
830119	LS/LS Press HT-2	5 x 3 g
830120	LS/LS Press HT-3	5 x 3 g
830121	LS/LS Press HT-4	5 x 3 g
	PORCELAIN 10 g–20 g	
830001	LS/Dentine D-A1	20 g
830002	LS/Dentine D-A2	20 g
830003	LS/Dentine D-A3	20 g
830004	LS/Dentine D-A3,5	20 g
830005	LS/Dentine D-A4	20 g
830006	LS/Dentine D-B1	20 g
830007	LS/Dentine D-B2	20 g
830008	LS/Dentine D-B3	20 g
830009	LS/Dentine D-B4	20 g
830010	LS/Dentine D-C1	20 g
830011	LS/Dentine D-C2	20 g
830012	LS/Dentine D-C3	20 g
830013	LS/Dentine D-C4	20 g
830014	LS/Dentine D-D2	20 g
830015	LS/Dentine D-D3	20 g
830016	LS/Dentine D-D4	20 g
830022	LS/Bleach Dentine BD-A	20 g
830023	LS/Bleach Dentine BD-B	20 g

Article no.	Article name	Content
830024	LS/Bleach Dentine BD-B0	20 g
830040	LS/Opaque Dentine OD-32	15 g
830041	LS/Opaque Dentine OD-37	15 g
830042	LS/Opaque Dentine OD-41	15 g
830043	LS/Opaque Dentine OD-43	15 g
830017	LS/Enamel E-57	20 g
830018	LS/Enamel E-58	20 g
830019	LS/Enamel E-59	20 g
830020	LS/Enamel E-60	20 g
830025	LS/Enamel Bleach S-AB	20 g
830021	LS/Clear CL-O	20 g
830026	LS/Transpa NT	20 g
 830056	LS/Transpa OT+	20 g
830057	LS/Transpa Geller-Mix	20 g
830027	LS/Transpa OT	20 g
830028	LS/Transpa TI-01	20 g
830029	LS/Transpa TI-02	20 g
830030	LS/Transpa TI-03	20 g
830031	LS/Transpa TI-04	20 g
830032	LS/Transpa TI-05	20 g
830033	LS/Pearl Enamel PS-0	20 g
830035	LS/Effect Enamel SI-02	20 g
830036	LS/Effect Enamel SI-04	20 g
830038	LS/Effect Enamel SO-10	20 g
830045	LS/Neck Transpa HT-51	15 g
830046	LS/Neck Transpa HT-52	15 g
830047	LS/Neck Transpa HT-53	15 g
830048	LS/Neck Transpa HT-54	15 g
830050	LS/Make In MI-61	15 g
830051	LS/Make In MI-63	15 g
830052	LS/Correction Powder KM	20 g
830055	LS/Glaze GL	10 g
LIQUIDS		
000016	Modelling Liquid low-fusing ML-L	25 ml
000017	Modelling Liquid low-fusing ML-L	50 ml
000018	Modelling Liquid low-fusing ML-L	250 ml
000010	Universal & Glaze Liquid UF/GL	25 ml
000011	Universal & Glaze Liquid UF/GL	50 ml
000090	Moisturing Fluid	25 ml

Article no.	Article name	Content
001210	Magic Colour Dentine Shade DS-A	3 g
001211	Magic Colour Dentine Shade DS-B	3 g
001212	Magic Colour Dentine Shade DS-C	3 g
001213	Magic Colour Dentine Shade DS-D	3 g
001214	Magic Colour Make In MI-61 ivory	3 g
001215	Magic Colour Make In MI-63 honey yellow	3 g
001216	Magic Colour Gingiva Shade GS-1 raspberry	3 g
001217	Magic Colour Gingiva Shade GS-2 salmon red	3 g
001218	Magic Colour Gingiva Shade GS-3 red violet	3 g
001221	Magic Colour Highlight HL-1 white	3 g
001222	Magic Colour Highlight HL-2 eggshell	3 g
001223	Magic Colour Highlight HL-3 lemon yellow	3 g
001224	Magic Colour Highlight HL-4 mandarine	3 g
001225	Magic Colour Highlight HL-5 flamingo	3 g
001226	Magic Colour Highlight HL-6 apricot	3 g
001227	Magic Colour Highlight HL-7 light brown	3 g
001228	Magic Colour Highlight HL-8 medium brown	3 g
001229	Magic Colour Highlight HL-9 olive	3 g
001230	Magic Colour Highlight HL-10 red	3 g
001231	Magic Colour Highlight HL-11 dove blue	3 g
001232	Magic Colour Highlight HL-12 grey	3 g
001233	Magic Colour Highlight HL-13 deep blue	3 g
001234	Magic Colour Highlight HL-14 fissure	3 g
001235	Magic Colour Highlight HL-15 illusion	3 g
001219	Magic Colour Glaze	3 g
001220	Magic Colour Glaze Fluo	3 g
001239	Magic Colour Glaze	10 g
001240	Magic Colour Glaze Fluo	10 g
	LIQUIDS	
001241	Magic Colour Liquid GL	25 ml
001242	Magic Colour Liquid GL	50 ml
001243	Magic Colour Special Liquid ML	25 ml
001244	Magic Colour Special Liquid ML	50 ml

Article no.	Article name	Content
001101	Make up Neo MN-F1 white	4 g
001102	Make up Neo MN-F2 eggshell white	4 g
001103	Make up Neo MN-F3 lemon yellow	4 g
001104	Make up Neo MN-F4 mandarine	4 g
001105	Make up Neo MN-F5 flamingo	4 g
001106	Make up Neo MN-F6 apricot	4 g
001107	Make up Neo MN-F7 light brown	4 g
001108	Make up Neo MN-F8 medium brown	4 g
001109	Make up Neo MN-F9 olive	4 g
001110	Make up Neo MN-F10 red	4 g
001111	Make up Neo MN-F11 blue	4 g
001112	Make up Neo MN-F12 grey	4 g
001113	Make up Neo MN-F13 deep blue	4 g
001114	Make up Neo MN-FA stain A	4 g
001115	Make up Neo MN-FB stain B	4 g
001116	Make up Neo MN-FC stain C	4 g
001117	Make up Neo MN-FD stain D	4 g
001121	Make up Neo Glaze GL-F	4 g
001151	In Nova Neo INN-1 melone	3 g
001152	In Nova Neo INN-2 rose	3 g
001153	In Nova Neo INN-3 yellow	3 g
001154	In Nova Neo INN-4 red brown	3 g
001155	In Nova Neo INN-5 brown	3 g
001156	In Nova Neo INN-6 olive	3 g
001157	In Nova Neo INN-7 grey	3 g
001158	In Nova Neo INN-8 blue	3 g
001159	In Nova Neo INN-9 white	3 g
001160	In Nova Neo INN-10 mais	3 g
001161	In Nova Neo INN-IL Illusion	3 g
001162	In Nova Neo INN-CRL Crackliner	3 g
001125	Make up Neo Fluid	25 ml
001126	Make up Neo Reconditioner Fluid	8 ml
001163	In Nova Neo Fluid	25 ml

CREATION LIQUIDS

Article no.	Article name	Content
000012	Modelling Liquid high-fusing ML-H	25 ml
000013	Modelling Liquid high-fusing ML-H	50 ml
000014	Modelling Liquid high-fusing ML-H	250 ml
000015	Modelling Liquid high-fusing ML-H	1000 ml
000016	Modelling Liquid low-fusing ML-L	25 ml
000017	Modelling Liquid low-fusing ML-L	50 ml
000018	Modelling Liquid low-fusing ML-L	250 ml
000019	Opaque Liquid high-fusing OL-H	25 ml
000020	Opaque Liquid high-fusing OL-H	50 ml
000021	Opaque Liquid high-fusing OL-H	250 ml
000010	Universal & Glaze Liquid UF/GL	25 ml
000011	Universal & Glaze Liquid UF/GL	50 ml
000025	Shoulder Liquid SPL	25 ml
000026	Shoulder Liquid SPL	50 ml
000090	Moisturing Fluid	25 ml
001241	Magic Colour Liquid GL	25 ml
001242	Magic Colour Liquid GL	50 ml
001243	Magic Colour Special Liquid ML	25 ml
001244	Magic Colour Special Liquid ML	50 ml
001125	Make up Neo Fluid	25 ml
001126	Make up Neo Reconditioner Fluid	8 ml
001163	In Nova Neo Fluid	25 ml



CREATION CERAMIC ACCESSORIES

Article no.	Article name	Content
000131	Brush Nr. I	1 pc.
000132	Brush Nr. II	1 pc.
000133	Brush Nr. III	1 pc.
000134	Brush Nr. IV	1 pc.
000135	Brush Nr. V	1 pc.
000136	Brush Kit	1 pc.
000321	Brush Creapast	1 pc.
000060	Sketch Pad	1 pc.
839912	LS/Manual DE-EN	1 pc.
910360	ZI-CT/Manual DE-EN	1 pc.
911020	CC/Manual DE-EN	1 pc.
911506	Ceramic meets Gingiva	1 pc.
911510	VC Brochure	1 pc.
911517	VC Workbook	1 pc.
911530	Magic Colour Manuel DE-EN	1 pc.
000059	CreaPen Separating Pen	1 pc.
000081	Creation Flüssige firing pad 10 g	1 pc.
000083	Creation Flüssige firing pad (3 x 10 g)	1 pack.

CREATION SHADE GUIDES

Article no.	Article name	Content
KIT-BASED SHADE GUIDES		
000064	CC/Shade Guide Kit 1 804	1 pc.
000068	CC/Shade Guide OM 11-17+G	1 pc.
000073	CC/Shade Guide Kit 2 1508	1 pc.
000087	CC/Shade Guide Transition Dentine	1 pc.
000092	ZI-CT/Shade Guide Kit 1	1 pc.
000093	ZI-CT/Shade Guide Kit 2	1 pc.
000076	ZI-F/Shade Guide Kit 1 (28 Shades)	1 pc.
000077	ZI-F/Shade Guide Kit 2 (26 Shades)	1 pc.
000058	Creation shade tab	1 pc.
000091	Shade Guide Wheel Empty	1 pc.
SHADE GUIDES		
000040	Shade Guide Dentine A-D (16 tabs)	1 pc.
000041	Shade Guide Opaque Dentine A-D (16 tabs)	1 pc.
000042	Shade Guide Transpa CL, UC, NT, OT (4 tabs)	1 pc.
000043	Shade Guide Opaque Dentine Modifier + AD OD-32, 37, 41, 43, 44, AD-1, 2 (7 tabs)	1 pc.
000044	Shade Guide Enamel E57-E60, SI01-SI06, SO10, SO11, PS0-PS3 (16 tabs)	1 pc.
000045	Shade Guide Transpa Modifier TI1-TI5, HT51-HT56 (11 tabs)	1 pc.
000046	Shade Guide Make In MI61-MI66 (6 tabs)	1 pc.
000047	Shade Guide Shoulder SP21-SP29 (9 tabs)	1 pc.
000048	Shade Guide Gingiva G1-G7, GN (8 tabs)	1 pc.
000049	Shade Guide Bleach BDA, BDB, BDB0, EAB, SP-AB (5 tabs)	1 pc.

Shades A1–D4 are referring to the Vita® Shade Guide.

Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

A Q U A L I N E

Article no.	Article name	Content
AMA-001	Aqualine master, mixing tray for ceramic, anthracite	1 pc.
AMA-002	Aqualine master, mixing tray for ceramic, white	1 pc.
AS-001	Aqualine smart, mixing tray for ceramic, anthracite	1 pc.
AS-002	Aqualine smart, mixing tray for ceramic, white	1 pc.
AM-001	Aqualine mini, mixing tray für Keramik, anthracite	1 pc.
AM-002	Aqualine mini, tray box with stain tray, with 15 cavities	1 pc.
AS-014	Aqualine Pool Set	1 pc.
AP-010	Aqualine master, glass plate, light grey	1 pc.
AP-011	Aqualine master, glass plate, dark grey	1 pc.
AS-010	Aqualine smart, glass plate, light grey	1 pc.
AS-011	Aqualine smart, glass plate, dark grey	1 pc.
AP-013	Aqualine master, cellulose strips (40 pcs.)	1 pack.
AS-012	Aqualine smart, cellulose strips (30 pcs.)	1 pack.
AM-010	Aqualine mini, cellulose strips (30 pcs.)	1 pack.
AMA-018	Aqualine master, labels, self-adhesive for stain (5 pcs.)	1 pack.
AMA-019	Aqualine master, labels, self-adhesive for ceramic (20 pcs.)	1 pack.
ST-003	Aqualine smart tray set, zircon stain plate (3 pcs.)	1 pc.
ST-004	Aqualine master tray set, zircon stain plate (4 pcs.)	1 pc.
AS-015	Aqualine Pool Set, pool set sponge (5 pcs.)	1 pack.
AC-001	Aqualine Glass Care (2 x 30 ml)	1 set
AC-002	Aqualine Two in One Build up Liquid (500 ml)	1 pc.

Notice



WILLI GELLER
Creation
SIMPLY BRILLIANT

YOUR LOCAL DEALER

Creation Willi Geller North America
3737 W 127th Street
0803 Alsip IL
USA

Phone +1 800 323 2164
info-northamerica@creation-willigeller.com
www.creation-willigeller.com/US

WILLI GELLER
Creation
SIMPLY BRILLIANT

DISTRIBUTOR
Creation Will Geller North America | 3737 W 127th Street, Alsip IL, 60803, United States
Phone +1 800 323 2164 | info-northamerica@creation-willigeller.com, www.creation-willigeller.com/US

MANUFACTURER
KLEMA Dentalprodukte GmbH | Koblacherstraße 3a, 6812 Meiningen, Austria